· Conscious dentistry ·

VIDU

EDTA GEL

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INSTRUCTIONS FOR USE

Description

A cleaning gel containing 17% Ethylene diamine tetraacetic acid. Used in canal irrigation and cleaning.

Composition

Ethylene diamine tetraacetic acid, poly ethylene glycol, distilled water, thickener agent

Indications

Used in In root canal treatment, in cleaned by washing of the root canal, in humidification of the root canal.

Soft and hard tissue chemically and physically removes residues.

Contraindications

Do not use EDTA GEL VIDU outside of the root canal system or on soft tissues. Immediately wash skin with soap and water after contact or flush eyes if contact occurs.

Application

After the isolation rubberdam entered into the root canal system. Pulp chamber is evacuated by an excavator and washed with irrigation solutions. Irrigation clinically applied with side disposable needles. Facilitated entry into the root canal using dental needle. The needle tip is placed in the root canal drained with moderate pressure EDTA GEL VIDU. Over flowing portion is impregnated with an aspirator. Canal file is placed in EDTA GEL VIDU full canal is continued cleaning and shaping process. Reapply if necessary until shaping complete. This process continues until the last use of the instrument. Keep off the syringe cover when not in use.

Side effects There are no known side effects

Storage

Protect from sunlight Store at 1-30 0C / 33,8-86 0F Keep in cool climates.

Cautions

Keep material out of children's reach!

For use in dentistry only!

Irritant be; Avoid contact with skin, eyes and mucous membranes. Disposal: Disposal must be made according to official regulations.

Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority.

Uncleaned packagings: Disposal must be made according to official regulations

Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes.

Expiry Date

From the moment the Syringe cover is opened should be used within 6 months. Manufacturing dates and LOT number on the top of the product packaging. Don't use after expiry date Applicator tips are single use.

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