· Conscious dentistry ·

DC FLOW CORE

en - Dual-cure material for core build-up and post cementation
es - Material de polimerización dual para la reconstrucción de muñones y la cementación de postes.
fr - Matériau à polymérisation double pour la reconstruction du moncone et la cementation di perrni
it - Materiale dual cure di reconstituzione coronaire e di scellament des tenons
pt - Material de polimerização dupla para enchimento com núcleo de cura dupla e cimentação para postes
de - Dual Material für Stumpfabauten und für die Zementierung von Wurzelkanalfittungen

INSTRUCTIONS FOR USE
PRODUCT DESCRIPTION: DC FLOW CORE is a dual curing radiopaque flowable microhybrid composite based material for core-build-up and cementation of posts. DC FLOW CORE is also suitable for cementing of crowns and bridges, inlays and onlays. For all indications the use of a dual cure bonding agent is required before application of the composite. DC FLOW CORE is based on poly- and difunctional methacrylates and inorganic filler particles of 0.05-2.5 µm. The total filler content is 66 % by weight and 46 % by volume. Delivered in auto-mixing 1:1 AUTOMIX-cartridges it can be easily dispensed and applied directly. DC FLOW CORE exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build-ups in cases where a light cure cannot be guaranteed to be sufficient. DC FLOW CORE meets the requirements of DIN EN ISO 4049, type 2, class 3.
INDICATIONS/INTENDED USE: Core build-ups, Cementing of posts, Cementing of crowns and bridges, inlays and onlays
PERFORMANCE FEATURES: The performance features of the product meet the requirements of the intended use.

CONTRAINDICATIONS: The placement of DC FLOW CORE is contraindicated if a dry working area or the recommended application technique are not possible. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

PATIENT TARGET GROUP: Persons who are treated during a dental procedure.
INTENDED USERS: This medical device should only be used by a professionally trained dental practitioner.

INCOMPATIBILITY WITH OTHER MATERIALS: Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

PREPARING OF THE AUTOMIX-CARTRIDGE: Take off the turn cap of the AUTOMIX-cartridge and throw it away (do not use it again). It is substituted by the special 1:1 mixing cannula by attaching it with a sideways turning movement, locking the cannula in place. The application gun is easy to load with the prepared cartridge and is ready for application.

Note: The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix. Store used cartridge with fixed used mixing cannula in the dark. The working time (23°C (74°F)) of DC FLOW CORE in the self cure mode is 1:30 minutes from start of mixing.

1. Post Cementation:
 - 1.1. Isolation:
Use of a rubber dam to isolate the tooth is strongly recommended.
 - 1.2. Root Canal Preparation
Refer to directions of the selected post manufacturer.
Before starting the preparation clean the tooth from residues. Prepare and clean the root canal with e.g. sodium hypochlorite solution, rinse and remove excess solution from the canal with a soft paper tip.
- Apply the bonding agent according to the corresponding instructions.
- 1.3. Post Cementation
Prepare the selected post according to manufacturer directions.

DC FLOW CORE is applied into the prepared root canal and onto the post. Seat the post carefully into the canal. Do not apply pressure until the post is seated. DC FLOW CORE self-cures within 3:30 minutes. For post stabilization light cure the coronal part of the cemented post for 20 seconds with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm2.

- As soon as the DC FLOW CORE has set proceed with the core-build-up procedure.
2. Core-Build-Up
 - 2.1 Isolation
Use of a rubber dam to isolate the tooth is strongly recommended.
 - 2.2. Cavity Preparation
Remove all existing old restorations and decay from the tooth. If necessary place any pins or posts. Refer to directions of the selected post manufacturer.

- 2.3. Pulp Protection
For pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.
- 2.4. Application of a Bonding Agent
Apply the bonding agent according to the corresponding instructions.
Note: It is essential that the primed dentine and enamel surfaces are dry and contaminant free for the application of DC FLOW CORE.
- 2.5. Application of DC FLOW CORE
Place the mixing cannula directly into the preparation and press out the paste.
DC FLOW CORE is automatically mixed when dispensed with slight and even pressure. Filling should occur from bottom upwards to prevent air voids. To facilitate placement of DC FLOW CORE place a matrix band around the prepared tooth.
DC FLOW CORE may be contoured by using a suitable modeling instrument.
Place DC FLOW CORE directly into the preparation and allow the system to self cure for 3:30 minutes. After that the material should be light-cured (40 seconds) with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm2. With this technique an optimum of physical properties will be obtained. An explnor can be used to test that DC FLOW CORE has completely set. Remove the matrix not earlier than the material has set.
Final core preparation on DC FLOW CORE can be carried out by using crown preparation burs.

ADDITIONAL NOTES/WARNINGS: The ambient light of the dental lamp may start

polymerization of the composite. Do not use any resin to adjust viscosity of the composite. Avoid contact with skin, mucous membrane and eyes.
If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required. Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates. Color stability meets the requirements of DIN EN ISO 4049. Commercial medical gloves do not protect against the sensitizing effect of methacrylates. Keep away from children!
COMPOSITION: Dimethacrylates, dental glass, silicone dioxide, photo initiators, catalysts
STORAGE: Do not store above 20 °C (68 °F). Protect from direct sunlight. Store unopened material in the refrigerator. Opened cartridges have to be used up within 3 months.
Do not use after expiry date.
DISPOSAL: Disposal of the product according to local authority regulations.
REPORTING OBLIGATION: Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.
Note: The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eu-damed>).
WARRANTY: The manufacturer warrants this product will be free from defects in material and manufacture. The manufacturer makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and the manufacturer's sole obligation shall be repair or replacement of the product.
LIMITATION OF LIABILITY: Except where prohibited by law, the manufacturer will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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