
Steam Sirio

- Instruction Manual
- Gebrauchs Anweisungen
- Manuale di Istruzioni
- Mode d'emploi
- Manual de manejo

CAUTION:

- **INSTALL THIS DEVICE IN A PROPER LOCATION**
- **TO BE USE BY AUTHORIZED AND INSTRUCTED PERSONNEL ONLY**

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as for directives:


**89/392/CEE
89/336/CEE
73/23/CEE
89/68/CEE**

we hereby certify under our responsibility that this device with matriculation:

nr. _____


is made in conformity with normatives:

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ELECTRIC DIAGRAM

PL – WIRE PLUG
 PD – FOOT PEDAL
 EV – ELECTRICAL STEAM VALVE
 LH – LAMP HOLDER WATER LEVEL
 HL – LAMP HEATING ELEMENT
 PS – PEDAL PANEL SWITCH
 GS - GENERAL SWITCH
 TH – THERMOSTATIC CONTROL
 PR - PRESSURE CONTROL
 HE – HEATING ELEMENT 1,5 KW
 CN – CONDENSATOR
 GR – GROUND
 FS – FUSE 10 AMP – T/10A 500V 6,3x32

Steam Cleaner

UNPACK

Verify integrity of the pack at the moment of the reception of the apparatus. Verify on the back side the yellow adhesive nameplate **QUALITY CONTROL PASSED** the control code is integral. Verify that in the packaging all the parts and the supplied accessories of series are contained together with the apparatus: steam gun with a set of 4 nozzles, loading funnel water, rod rest-gun, electrical pedal, user's documentation and warranty.

INSTALLATION

Take the unit out of the packaging and place it on a flat surface strong enough to bear its weight (14 Kg. with water). Connect the electrical pedal (ref. 2) to the appropriate outlet placed on the upper side of the unit. Insert the rod rest-gun and screw it in the appropriate site (ref. 40). Connect the steam tube (ref. 1) to the appropriate sleeve of electric valve (ref. 9) and properly fix it with the equipped band. Connect the main plug to a power outlet.

ATTENTION: make sure that the voltage corresponds to the one brought back on the fabrication nameplate placed on the back side of the apparatus. Press the power switch on the left to position 1 (ref. 31 SX) to turn on the unit and the inside water level light. Remove the filling stopper (ref. 8) making a light pressure downwards; then fill the boiler with demineralized water, untill you reach the maximum level showed on the front panel of the unit. Screw the stopper (ref. 8) on the hole of threaded loading to close. Now turn off the left switch and turn on the switch on the right (ref. 31 DX) to start the heating; the ignition of the heating element is marked by the red light (ref. 33) placed on the left side of both the two switches. When the pressure of 3.5 bars is reached the red light automatically switch off, now the unit is ready to work. Possible boiling noise when starting can be eliminated letting out the inside air bubbles by pressing the steam spilling pedal. To obtain the steam out of the gun (ref. 1), turn on the switch on the left (the same which turns on the water level light), then press the electrical pedal (ref. 2) previously connected. Regulate the steam flow by rotating the regulation grip handle placed on the base of the electric valve (ref. 9). The initial steam capacity could be very short, wait for a few minutes till the pressure is reached again, then the steam will be released continuously. When you end steaming, turn off the switch on the left to disconnect the pedal and electric valve. At the end of the day turn off the main switch on the right, without unloading the pressure.

MAINTENANCE

Check the level of the water: verify every time the unit is turned on that the level of the water is not under the minimum level reference showed on the front panel. If the level of the water is too low, fill the boiler following the procedure of loading of the water (see points 6-7-8).

- **Drainage of the residual calcium in the boiler** - Every 6 months of operating, the drainage of the residual calcium in the boiler must be carried out executing the following operations: control that the apparatus is cold and without internal pressure. Disconnect the power cord from the electrical outlet (rif. 39). Place the apparatus on a bench on its left side. Unscrew the nut placed under the base of the boiler (rif.54) without removing it completely to avoid to bathe the inside of the apparatus. Place under the apparatus a washbasin or a container of at least 5 liters capacity. Then unscrew completely the nut from the boiler and remove it from the center. Remove now the stopper (rif. 8) and make discharge the residual water completely. Pour inside with the funnel 3/4

liters of water for rinsing the boiler. Then screw again not too strongly the nut under the boiler, placing a tape of teflon on the thread (approximately 8/9 complete turns). Reconnect now the apparatus to the power cord, and fill the boiler with demineralised water as usual.

GUIDE TO COMMON QUESTIONS

The apparatus is not turning on:

1. Control that the power cord plug is properly connected to the electrical outlet.
2. Check that power is on at the electrical outlet.

Steam flows little in relation to the pressure indicated on the gauge:

Check that the regulation grip handle of the steam placed on the electric valve (rif. 9) is completely opened.

A strong water boiling noise from the boiler when start heating (rare operation):

Open slowly the filling cap to let air bubbles out.

After reached the pressure the steam comes out for limited time:

The machine's components are cold, please allow necessary time to warm up all internal components (15-20 minutes at 3.5 bar without release of steam)

Carried out these controls, in case the apparatus still presents disadvantages, contact the customers' technical service.

SPECIFICATIONS

W/Hz 230/50 - 60 or 110/50

Absorption 1500 Watt

Weight Kg. 10

Internal Fuse 10 AMP – In case of substitution of the fuse ask for technical assistance.

DIMENSIONS

Height cm. 50 Width cm. 22 Depth cm. 30

NOTES

For the external cleaning of the device use a cotton tawel damped with water.

The device can work in laboratory room from 0°C to 35°C.

Noise level (1 mt): stand-by 26.6 db, warm-up 37.5 db, steam cleaning 68.7 db.

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