R-100247

E.....	HORNO DE CERAMICA	2
GB	CERAMIC FURNACE.....	11

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

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NOTE: Carefully read this manual. It contains important safety information in the installation, use and maintenance. If the equipment is used in a manner not specified by the manufacturer, the security provided by the equipment may be compromised.

Meaning of the symbols on the unit

Danger. High temperature.


Unpacking

- Remove the parking belt and carton.
- Remove the plastic bag.
- Remove the outsider foam.
- Take out the firing platform and firing plate which are placed in the upper foam.
- Carefully take out the furnace from foam base. Place it on a stable surface.
- Take the vacuum pump out of the carton. Place it on the ground nearby the furnace.
- Keep the original packaging for future transport purposes.

Vacuum pump installation

- Connect the vacuum pump hose to the vacuum hose connector at the rear of the furnace.
- Plug the vacuum pump power line into vacuum pump electrical outlet at the rear of the furnace.

1. Solenoid valve.
2. Filters.
3. Exhaust hole.
4. Hose.
5. Power cord.

Furnace installation

- Connect the plug to a 230 V, 50/60 Hz earthed power outlet.
- Press power switch ON at the rear of the furnace. It will enter the stand by menu.
- Press ↑ to lift the muffle.
- Place the firing platform on the middle of the tray. Make sure that it is placed steadily.
- **After switching the furnace on, the firing platform must be immediately placed on the tray; otherwise the sealing ring on the tray will be burn out. The firing platform can not be taken out from the furnace except for transport.**

Precautions

- Keep these instructions for further consultation.
- Check out the state of the unit upon unpacking it.
- Damage caused by a faulty installation or a wrong use of the unit is not covered by the warranty.
- To prevent electric shock hazards, connect to a proper power source (230 V, 50/60 Hz, earthed). Before connecting the unit, make sure the voltage shown in the identification plate matches that of the power source to be used. The installation must fulfil all the applicable safety standards and the instructions of the manufacturer. The manufacturer will not be responsible of the damage caused by an incorrect ground connection.
- If the power supply is not stable, the furnace will be unable to work normally.
- If sudden power failure when the furnace is working, unplug it.
- Remove any dirt or water remains from the cord connection.
- Do not pull, twist or fold the power cord excessively.
- Do not touch the unit or the power cord with wet hands.
- Do not install the unit near to heat sources, dust, high humidity or near inflammable materials (alcohol, solvents,...).
- Unplug the unit immediately whenever you notice an unusual noise, smell or smoke.
- Do not allow children or non-qualified personnel handling the unit.
- When in doubt, do not start the unit. Contact the supplier for questions or detailed information about problems with the unit.
- Do not dismount, repair or modify the machine.
- Unplug the unit in case of storm and when you are not going to use it for a long time.
- Always the furnace is used, the surface of the furnace and the muffle should keep clean.
- Clean the furnace only with a dry or slightly damp cloth.


- Unplug the unit before cleaning, do not. Do not wet the unit. Do not use chemical products to clean it.
- Do not place sundry objects around the furnace.
- Before using the furnace, make sure it is on a firm and strong surface.
- Before switching the furnace off, press ↓ to make the muffle move down.
- **The moisture in the muffle should be correctly removed, otherwise the moisture may enter the vacuum system to damage the vacuum pump and solenoid valve etc. Press key [PURGE].**

Intended use

The R-100247 is a ceramic furnace especially designed for the field of dental mechanics

Packing contents

- 1 vacuum pump.
- 1 firing platform.
- 1 calibration kit.
- 1 tray.


1. **Display screen**
2. **Numeric keys (0-9)**
3. **CANCEL:** When editing parameters, one figure of parameter value will be cancelled.
4. **ENTER:** To confirm parameter and return to stand-by menu.
5. **PURGE:** To activate the recommended program if the muffle is moistened, contaminated or not used for a long time.
6. **SET:** To enter into the System Setting Menu.
7. **STOP:** To stop all current programs including the muffle movement, firing programs, etc.
8. **START:** To start the selected firing program.
9. **← (left) and → (right):** To select the parameters of the firing program.
10. **↑ (up) and ↓ (down):** To move the muffle. To select parameters in System Setting Menu.

Menú stand-by


After switching the furnace on, the screen will show the stand-by menu.

1. Date and time.
 2. Program number.
 3. Firing curve and parameter values.
 4. Room temperature.
 5. Muffle temperature.
- In stand-by menu, press ← and → to select the different parameters.
 - Each time you press ← and →, a parameter will highlight. At the same time, the parameter name will be displayed on left lower corner.
 - In stand-by menu, press ↑ and ↓ to move the muffle up and down.
 - In stand-by menu, when **START** is pressed, the firing program begins to run and the screen will change to the RUN menu.
 - In stand-by menu, when numeric keys (0-9) are pressed, the firing program will be selected (P0001-P200). Press **ENTER** to confirm the program number and return to stand-by menu.


Select firing program number

- In stand-by menu, directly press numeric keys (0-9) to select a program number.
- The program number (P001-P200) will highlight.
- Press **ENTER** to confirm program and return to stand-by menu.


→


Editing firing cycle parameters

- In stand-by menu, press ← and → to select the different parameters.
- Each time you press ← and →, a parameter will highlight. At the same time, the parameter name will be displayed on left lower corner.
- When a parameter is highlighted, press numeric keys (0-9) to edit the value.
- Press **ENTER** to save the parameter. The following parameter will highlight.
- After setting parameters, when **START** is pressed, the parameter is automatically saved. The firing program will begin to run and the screen will change to the RUN menu.
- The unit of time is second, the unit of temperature is degree Centigrade (°C).
- If the vacuum level is more than 97kPa, the vacuum pump will work automatically until vacuum is released.


Idle (100 - 750 °C)

The highest temperature that the muffle will reach and keep in stan-by mode.


Dry time (0 - 600 seconds)

The length of time that the muffle moves to the lower position alter the muffle temperature reaches to Idle temperature.


Preheat time (0 - 600 seconds)

The length of time that the muffle moves to the bottom position after dry time is finished.


Heat rate (2 – 140 °C/min.)

The rate at which muffle temperature increases with time.


Hi temp (410 - 1305 °C)

The maximum temperature during the firing program
NOTE: If *Hi Temp* is lower than *Idle* temperature, the furnace will automatically revise *Hi temp* to be 10 °C higher than *Idle* temperature.


Hi Temp. Hold (0 - 900 seconds)

The length of time that the temperature will be held after the muffle temperature reaches to Hi temp.
NOTE: in order to ensure that chamber can not open until the vacuum is released in the muffle, if *Hi Temp. hld* time is less than *Vac. Hold* time, the furnace should revise *Vac. Hold* time to 30 seconds shorter than *Hi Temp. Hold* time.


Cool time (0 - 600 seconds)

The length of time that the muffle will keep down after Hi Temp. hold.

Vac. level (0 - 100 kPa)

The value of vacuum which will be maintained in the muffle.


Vac. start

The temperature at which vacuum pump will begin to run.

Minimum: \geq Idle temperature.

Maximum: 10 °C less than HI Temp


Vac. stop

The temperature at which vacuum in the muffle begins to release.

Minimum: 10 °C more than Idle temperature.

Maximum: \leq Hi temperature.

Vac. hold

The length of time that vacuum in the muffle will be maintained after the muffle temperature reaches to Hi Temp.


It must be set to 30seconds less than Hi Temp Hold time.

Vacuum

When Vac. Hold time is set to 0, Vac. Start and Vac. Stop must be edited.

For example:


When VAC. HOLD is 0 second, if VAC. STOP temperature is 930 °C, the vacuum in the muffle will be released after the muffle temperature reaches 930 °C.

When Vac. Hold time is set to > 0 seconds.


For example: in petal-porcelain program.

If VAC. HOLD time is edited to 30 seconds, when muffle temperature reaches to HI TEMP., the vacuum in the muffle will be released after 30 seconds.


Menú RUN

- In Stand-by Menu, if RUN is pressed, the firing program will begin to run. At the same time, the screen will change to RUN Menu. The progress of the firing cycle program is shown by a curve.
- While the firing cycle runs, the actual step will highlight.
- In RUN menu, if STOP is pressed, the firing program will finish.
- After the firing cycle is completed, the screen will return to stand-by menu.
 1. Date and time.
 2. Program number.
 3. Firing curve and parameter values.
 4. Room temperature.
 5. Muffle temperature.

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
Date and time


- Press ↑ and ↓ to select the parameter Time. The menu of time and date will be shown.
- Press ← and → to select year, month, day, hours, minutes and seconds. The parameter will highlight.
- Edit each parameter by numeric keys (0-9).
- After setting the date and time, press ENTER to save and press ENTER again to return to stand-by menu.

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
- Press ←. Display will show symbol – (decrease the muffle temperature).
- Press →. Display will show symbol + (increase the muffle temperature).
- Use the numeric keys (0-9) to edit the temperature value that you want to increase or decrease on the basis of original temperature.
- Then press ENTER to save and press ENTER again to return to stand-by menu.


Language

- Press ↑ and ↓ to select the parameter Language.
- Press ← and → to choose the language.
- After setting, press ENTER to save and press ENTER again to return to stand-by menu.

Changing vacuum unit

- Press ↑ and ↓ to select the parameter Vac. unit.
- Press ← and → to choose kPa, mmHg or inHg.
- After setting press ENTER to save and press ENTER again to return to stand-by menu.

Introduction of programs

Programs 001–010

- Programs 001–010 are without vacuum protection. They can be used if vacuum is not need or not required strictly (for example when vacuum system has trouble).
- When leaking in the vacuum system after setting vacuum level is ready, the vacuum pump will not begin to run automatically.
- When the vacuum level in the muffle is less than SETTING VACUUM, the furnace will also not stop working.

Programs 011–100

- Program 011–100 are the programs with partial vacuum protection.
- They can be used if a high vacuum level is need at the initial stage of firing program but not required strictly afterwards (for example when the sealing of the muffle is not well).
- When the vacuum level in the muffle doesn't reach SETTING VACUUM, the warning tone will sound and the furnace will stop working.
- When the vacuum level in the muffle reaches to SETTING VACUUM, if the vacuum level in the muffle is decreased because of leaking, the vacuum pump will not start working again and the furnace will also not stop working.

Programs 101–200

- Program 101–200 are the programs with full vacuum protection.
- When the vacuum level in the muffle doesn't reach to SETTING VACUUM, the warning tone will sound and the furnace program will stop working. It needs to switch off the furnace and check it. The furnace can't use again until troubleshooting.
- After the vacuum level in the muffle reaches SETTING VACUUM, if the vacuum level is decreased because of leaking, the vacuum pump will start working again to make the vacuum level reach to the SETTING VACUUM again.
- Program No.121 and 122 is pre-set the Ivoclar pressing program. It is only for reference. When Program 121 begins to run, it will automatically enter Program 122 until finish whole program. Program No.122 can't be used separately.

Alarm messages

Thermocouple alarm

The message Thermocouple Fault will appear if:

- Thermocouple is broken.
- Thermocouple is disconnected.
- Temperature measuring circuit is faulty.

Heather alarm

The message Heater Fault will appear if:

- Heater coil burns out.
- Lead connected with Heather coil is broken.
- Heating circuit is faulty.

Vacuum alarm

The message **Vacuum Fault** will appear is vacuum fails. Check vacuum and seal system.

Limit switch alarm

The message **Limit switch Fault** will appear if the limit switch is damaged. Check limit switch inside the furnace.

Temperature calibration


If user finds that muffle temperature is on the high or low side Turing operating furnace, it can be readjusted within the range of 0 to 80 °C.

NOTE: the temperature calibration test should be used alter the muffle is completely dried. Negligence will cause incorrect temperature calibration.

- Cut a 1-2 mm length silver sample and put it in the calibration set crucible. Then place the crucible on the middle of the firing plate.
- Select program n. 8 with the following parameters:

Dry time:	0 s
Preheat time:	0 s
Firing temperature:	955 °C
Firing time:	180 s
Cool time:	0 s
Vacuum:	0 kPa
Heat rate:	60 °C/min
- After setting parameters, start this program.
- After this program runs, the silver sample should not have been fused. If it fused, it means that the muffle temperature is too high. Calibrate the temperature as previously indicated and repeat the test until the sample does not fuse.
- One the sample does not fuse, select program n. 9 with the following parameters:

Dry time:	0 s
Preheat time:	0 s
Firing temperature:	965 °C
Firing time:	180 s
Cool time:	0 s
Vacuum:	0 kPa
Heat rate:	60 °C/min
- After setting parameters, start this program.
- After this program runs, the silver sample should have been fused. If it did not fused, it means that the muffle temperature is too low. Calibrate the temperature as previously indicated and repeat the test until the sample fuses.

Original shape of the sample before calibration.

Only the middle part of the sample is bulged.

Silver sample is melted to the round.

Vacuum conversion table

inHg	kPa	mmHg	inHg	kPa	mmHg
1	4	20	16	53	400
2	7	50	17	57	430
3	10	70	18	60	450
4	13	100	19	64	480
5	17	120	20	68	510
6	21	150	21	71	530
7	24	170	22	74	560
8	27	200	23	78	580
9	30	220	24	81	610
10	34	250	25	85	630
11	37	270	26	88	660
12	40	300	27	91	680
13	43	320	28	94	710
14	47	350	29	98	740
15	50	370	30	101	760

1 inHg = 3.38 kPa = 25.4 mmHg

1 MPa = 1000 kPa = 1000000 Pa

Technical characteristics

Electrical supply: 230 ± 22 V, 50/60 Hz, 15 A
 Power: < 1700 W
 Height: 380 mm
 Width: 340 mm
 Depth: 435 mm
 Chamber: Ø9 x 6 cm
 Weight: 26 kg (no pump)
 Maximum temperature: 1200 °C
 Maximum vacuum: -96 kPa