· Conscious dentistry ·

UNIVERSAL COMPOSITE

- en** · Light-curing universal hybrid composite
- es** · Composite híbrido universal fotopolimerizable
- it** · Composito universale ibrido fotopolimerizzabile
- fr** · Composite hybride universal photopolymérisable
- pt** · Compósito híbrido universal fotopolimerizável
- de** · Lichthärtendes universal Hybrid Composit

"THE FUTURE
IS IN YOUR HANDS"

VIDU UNIVERSAL COMPOSITE

en

Instructions for use

Patient target group
Children, adolescents, adults and older patients

Intended users
Dentist

Composition
Qualitative data: oligomeric methacrylates, silicon dioxide, glass ceramics
Quantitative data in % G/G: Oligomeric methacrylates 23, fillers mid 4, mini 65 and nano 8

Intended purpose of the product
UNIVERSAL COMPOSITE VIDU is a light-curing, radiopaque universal filling composite.

Indications
Restoration of all cavity classes.

Contraindications
Use of UNIVERSAL COMPOSITE VIDU is contraindicated if drainage or the recommended application method is not possible. Do not use directly on opened pulp or if there is a known allergy to any ingredient.

Performance characteristics of the product
UNIVERSAL COMPOSITE VIDU is a universal filling composite that cures with light.

Specific precautions
None

Application
1. Clean the tooth before preparation, remove all cleaning residues with water spray.
Determine the Vita® tooth colour. Using a dental dam is recommended.
2. Prepare the cavity according to the general rules of adhesive technique. Bevel the edges of the enamel and protect areas close to the pulp.
3. If required, etch according to your preferred technique. Observe the instructions for the materials used.
4. Apply adhesion promoter
5. If necessary, place an underfill of thin flowing composite or Compomer. Using a die is recommended.
6. Fill the cavity in the desired colour with layers no thicker than 2 mm.
Leave each layer to light cure for 20 seconds.
7. Remove excess, contour and finish e.g. with diamond finishers and flexible discs of decreasing grain size. A high-gloss finish can be achieved with polishing pastes and polishing pads/discs.

Processing instructions
The processing instructions and precautionary measures described in the instructions for use accompanying the product must always be observed. If, despite adherence to these instructions, an adverse reaction occurs in the patient or other adverse events occur, we request that you inform us in any case, describing the accompanying circumstances as precisely as possible.

Warnings
Residual risks
Please observe carefully
• Light units should be emitting at 450 nm and should be checked regularly. Place the light as close as possible to the filling material.
• Do not remove the oxygen-inhibited layer, as this is needed for bonding with the next layer.
• Follow the instructions for use of other materials used with UNIVERSAL COMPOSITE VIDU.

NOTE
The application aids included with the product (cannulas, disposable brushes, etc.) are only for single use on a single patient!

Side effects
No side effects have been reported to date. In individual cases, hypersensitivity to components of the material may occur.
Avoid skin contact with the paste. In case of accidental contact, wash immediately with water and soap. In case of eye contact, rinse immediately with plenty of water and consult a doctor. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible.

Risks in case of overdose are not known.

Interactions
Avoid materials containing eugenol, moisture and oily air, as they can prevent polymerisation at the point of contact. Discolouration may occur when using cationic mouthwashes, plaque disclosers or chlorhexidine.

Storage
The product should be stored between 0°C - 22°C. After use, empty the syringes immediately so that the material is not hardened by the ambient light.

Minimum shelf life
Do not use after the expiry date

Variants
REF VI92001: UNIVERSAL COMPOSITE, refill, A1, 4,5 g
REF VI92002: UNIVERSAL COMPOSITE, refill, A2, 4,5 g
REF VI92003: UNIVERSAL COMPOSITE, refill, A3, 4,5 g
REF VI92004: UNIVERSAL COMPOSITE, refill, A3.5, 4,5 g
REF VI92005: UNIVERSAL COMPOSITE, refill, B1, 4,5 g
REF VI92006: UNIVERSAL COMPOSITE, refill, B2, 4,5 g
REF VI92008: UNIVERSAL COMPOSITE, refill, C2, 4,5 g
REF VI92010: UNIVERSAL COMPOSITE, refill, D2, 4,5 g
REF VI92022: UNIVERSAL COMPOSITE, refill, A2 opaque, 4,5 g
REF VI92023: UNIVERSAL COMPOSITE, refill, A3 opaque, 4,5 g
REF VI92019: UNIVERSAL COMPOSITE, refill, Incisal, 4,5 g
REF VI92020: UNIVERSAL COMPOSITE, refill, Bleach, 4,5 g

Disposal
Unused material can be disposed of with household waste after curing.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

  H317 May cause allergic skin reactions
H318 Causes serious eye damage
H411 toxic to aquatic organisms, with long-lasting effect
Danger **Danger**

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