ULTRASONIC SCALER TIP

en • Ultrasonic Scaler Tip **es** • Punta de Ultrasonido **fr** • Punta per Ultrasouni **it** • Inserti à Ultrasoni **pt** • Ponta Ultrassonica **de** • Ultraschall Spitze

INSTRUCTIONS FOR USE en

Product Name: Ultrasonic Scaler Tip **Model:** (See label on packaging) **Product Description:** The worktip is replaceable and attaches to the handpiece of the ultrasonic scaler. Fed by the ultrasonic scaler, the tip or work insert is usually made up of the workpiece and the shank. **Intended use:** The tip of the ultrasonic scaler coupled to the ultrasonic scaler handpiece is used to clean and shape the tooth surface and root canal. **Patient Population and Indications:** Adults and pediatrics with the following indications: Regular periodontal cleaning. Treatment of periodontal disease. **Intended Users:** The ultrasonic scaler is designed to be used by qualified and professionally trained dentists.

Contraindications:

- Patients with hemophilia are prohibited from using this equipment.
- Patients or doctors with cardiac pacemakers are prohibited from using this equipment.
- Patients with heart disease, pregnant women, and children should exercise caution when using the equipment.

Compatibility Descriptions:

- Gx, GxT (x:1-999) Female Scratch Screw Series: Compatible with EMS Brand 1. Px, PxL, PxR, PxLT, PxRT, PxLDT, PxRDT, PxLD, PxRD, PxD, PxT, PxDT (x:1-999) Female screw series for perio: compatible with EMS brand ultrasonic scalers.
- Ex, ExD, EXT, ExDT (x:1-999) Female Screw Series for Endo: Compatible with EMS brand ultrasonic scaler.
- GDx, GDxT (x:1-999) Female Scraping Screw Series: Compatible with SATELEC and NSK ultrasonic scaler.
- PdX, PDXD, PDXT, PDXL, PDXR, PDXLT, PDXRT, PDXLDT, PDXRDT, PDXLD, PDXRD, PDXDT (x:1-999) Female screw series for Perio: compatible with SATELEC and NSK ultrasonic scaler.

- EDx, EDxD, EDXT, EDxDT (x:1-999) Female Screw Series for End
- GSx (x:1-999) Male Scraping Screw Series: Compatible with Sirona ultrasonic scaler.

- PSx, PSxD (x:1-999) Male screw series for Perio: compatible with Sirona ultrasonic scalers.
- ESx, ESxT, ESxD, ESxDt (x:1-999) Male screw series for Endo: compatible with Sirona ultrasonic scalers.
- GCx, GKx (x:1-999) Male Scraping Screw Series: Compatible with KAVO Climber.

- EKx (x:1-999) Male Screw Series for Endo: Compatible with KAVO Climber.
- Ax, AxT (x:1-999) Amdent series scraper tip, compatible with Amdent ultrasonic scrapers.
- SBx, SBxT, SBxR, SBxL, SBxRT, SBxLT (x:1-999) Female Screw Series for Cavity Preparation: Compatible with EMS ultrasonic scaler.
- SBDX, SBDxT, SBDxR, SBDxL, SBDxRT, SBDxLT (x:1-999) Female screw series for cavity preparation: compatible with SATELEC and NSK ultrasonic scarifiers.
- SBSx(x:1-999),SBSL, SBSR: Male Screw Series for Cavity Preparation: Compatible with SIRONA Ultrasonic Scalers.
- USx, USxL, USxR, ULx, UCx, UPx, UEx, U(x:1-999) Oral surgery series tips, compatible with Mectron bone surgery devices.
- 201, 202, 203: Female scraping screw series, compatible with KAVO ultrasonic scrapers.
- ETx (x:1-999) Female screw series for Endo, compatible with Satelec ultrasonic scaler.
- ExD-S (x:1-999) Female screw series for Endo, compatible with NSK ultrasonic scaler.
- Sx (x:1-999) Female scraping screw series, compatible with NSK and Amdent ultrasonic scaler.
- PCx(x:1-999) Female Scraping Screw Series, Compatible with KAVO Ultrasonic Scaler.

How to use:

- Scraping and cavity preparation tips should be used with the scraping function of the ultrasonic scraper. Perio tips should be used with the climber's Perio function. Endo tips should be used with the climber's Endo function.
- Select the appropriate tip according to the operating requirements, firmly install it on the handpiece with a torque wrench.
- Turn on the ultrasonic scaler for operation. During the operation, do not allow the tip of the appliance to come into contact with the surface of the teeth in a vertical direction, do not exert too much pressure so as not to damage the teeth or the tip.

- After operation, rinse the tip and handpiece, leaving the appliance running for 30 seconds with water, then remove the tip and handpiece for sterilization.

Precautions:

- Make sure the tip is tightly fitted to the handpiece and spray well during operation.
- The climber's tip should be cleaned, disinfected, and sterilized before each use. Follow the recommendations in the manual in the appendix "Cleaning, Disinfection and Sterilization Reprocessing Instructions".

The following sterilization methods are prohibited:

- Balls in water.
- Immersion in iodine, alcohol and glutaraldehyde.
- Place in the oven or microwave.
- The maximum cleaning, disinfection and sterilization cycles for the tip of the ultrasonic scraper are identified as 300 cycles.
- Do not bend or polish the tip.

- After the gold-coated tips have come into contact with an oxidizing solution such as H2O2 or NaClO, use clean water to rinse the tip.

Storage: The tip should be stored in a clean, dry and ventilated indoor environment with a relative humidity of 10% to 93%, an atmospheric pressure of 70kPa - 106kPa, a temperature of -20°C + 40°C without corrosive gas.

Environmental protection: The product does not contain batteries or toxic substances. And there are no components that need to be specially removed from the main unit for disposal. Once the life of the appliance is over, do not dispose of it in the household waste. Comply with your country's Waste Electrical and Electronic Equipment (WEEE) directives and medical waste disposal regulations. The tips, endochuck, torque wrench, endocanal wrench and handpiece, which can easily come into contact with biological sources and cause biological hazards, must be separated from the main unit and reprocessed before being discarded and destroyed. Tips and endochucks are sharp instruments that can easily scratch people. They should be disposed of in hospital waste bins for sharp instruments. Once the life of the device is over, do not dispose of it in the household waste. Comply with your country's Waste Electrical and Electronic Equipment (WEEE) directives and medical waste disposal regulations. The tips, endochuck, torque wrench, endocanal wrench, and handpiece, which can easily come into contact with biological sources and cause biohazards, must be separated from the main unit and reprocessed before being discarded and destroyed. Tips and endochucks are sharp instruments that can easily scratch people. They should be disposed of in hospital waste containers for sharp instruments.

Cleaning, Disinfecting, and Sterilizing Reprocessing Instructions:

- Start of Work
- Please read these instructions carefully as they explain all the most important details and procedures. Pay special attention to safety precautions. Keep these instructions handy at all times.

- In order to avoid personal and material damage, it is necessary to take into account the relevant Directives.
- The instructions contained in this manual apply only to the product with which it was delivered.

- Introduction
- 1 The these reprocessing instructions provide instructions for cleaning, disinfecting, sterilizing, and packaging the manufacturer's reusable products intended for reprocessing in medical facilities.
- 2 The objective of reprocessing reusable products is to reduce the bioburden and achieve sterility of these products in order to eliminate the risk of infection related to the reuse of the product. Decisions regarding the cleaning, disinfection or sterilization of the manufacturer's medical and dental instruments are based on the potential risk of infection associated with their use.
- 2.3 The use of steam sterilization is recommended.
- 2.4 Remember that a high level of sterilization or disinfection cannot be achieved without first cleaning the elements of the assembly.
- 2.5 If you feel that the manufacturer's reprocessing instructions appear to be inadequate, please inform the manufacturer of these inadequacies.
- 2.6 You are encouraged to report adverse events related to the reprocessing of products. Report these developments directly to the manufacturer.

- Reprocessing: Instructions for Reusable Products
- 3.1 The instructions are binding for the reprocessing of all reusable products (hereinafter referred to as "products") of the manufacturer. If necessary, additional product-specific instructions are included.

Important: Before use, carefully read the instruction manual of the manufacturer of the instrument and the devices with which the product will be used.

- 3.2 Reusable products must be cleaned, disinfected and sterilized before first use. Reprocessing procedures have only limited implications for this device. Therefore, the limitation of the number of reprocessing procedures is determined by the function/wear of the product. In terms of processing, a maximum number of reprocesses is not allowed. The device should no longer be reused if there are signs of material degradation.

In case of damage, the product must be reprocessed before being sent to the manufacturer for repair.

4. Preparation - basic principles
- 4.1 Effective sterilization can only be carried out after effective cleaning and disinfection has been completed. As part of your responsibility for the sterility of products during use, ensure that only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that validated parameters are adhered to during each cycle.

- 4.2 Also respect the legal requirements applicable in your country, as well as the hygiene standards of the hospital or clinic. This applies especially to the additional requirements for prion inactivation.
5. Preparation at the point of use: Turn off the product. Remove coarse dirt from products with cold water (<40°C) immediately after use. Do not use setting detergent or hot water (>40° 0), as residues may be deposited that may influence the outcome of the reprocessing process.

- Store products in a humid environment
6. Transportation: Store and transport safely to the reprocessing area to avoid any damage and contamination of the environment.
7. Preparing for Decontamination: Products should be reprocessed, as far as possible, disassembled.
8. Pre-cleaning: Pre-clean products manually until they are visually clean. Soak the products in a cleaning solution and rinse the lumens with a cold tap water jet gun for at least 10 seconds. Clean surfaces with a soft-bristled brush.
9. Cleaning: With regard to cleaning/disinfection, rinsing and drying, a distinction should be made between manual and automated reprocessing methods. Preference should be given to automated reprocessing methods, especially because of their greater potential for standardisation and industrial safety. Automated Cleaning: Use a washer disinfector (WD) that meets the requirements of the ISO15883 series. Place the products in the machine on a tray. Connect the products to the WD using the appropriate adapter and start the program:
 - 4 min pre-wash with cold water (<40°C)
 - Emptying
 - 5 min minimum rinse with a mild alkaline cleaner at 55°C
 - Emptying
 - 3 min neutralization with hot water (>40°C)
 - Emptying
 - 5 min intermediate rinsing with lukewarm water (>40°C)
 - Emptying
- Automated cleaning processes were validated using 0.5 neodisher MediClean forte (Dr Weigert).

According to EN ISO 17664, manual reprocessing methods are not required for these devices. If a manual reprocessing method must be used, validate it before use.

10. Disinfection: Automated thermal disinfection in the washing machine/disinfector, taking into account national requirements regarding the A0 value (see EN 15883). A 5-minute disinfection cycle at 93°C has been validated to achieve an A0 value of 3000.

11. Drying: Automatic drying:
 - External drying of the products at 40° C, 5 min using the drying cycle of the washer/disinfector. If necessary, additional manual drying can be done using a lint-free towel. Inflate product cavities using sterilized compressed air.
12. Functional Testing, Maintenance: Visual inspection of product cleanliness and reassembly if necessary. Functional test according to the instructions for use. If necessary, rewash again until the products are visibly clean. Before autoclaving and sterilizing, make sure that the products have been kept in accordance with the manufacturer's instructions.

13. Packaging: Pack the products in packaging material suitable for sterilization. The material and packaging system conform to the EN ISO11607 standard. The manufacturer assumes no responsibility for the use of other sterilization procedures (e.g., ethylene oxide, formaldehyde, and low-temperature plasma sterilization). In these cases, it is necessary to observe the respective valid standards (EN ISO 14937/ANSI AAMI ISO 14937 or the specific standard for the procedure) and check the suitability and effectiveness of the procedure in principle (if necessary, including investigations on residues of sterilising agents), taking into account the specific geometry of the product as part of the validation.

- Maximum sterilization temperature 138°C.
14. Sterilization: Sterilization of products by applying a fractional pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665), taking into account the requirements of the respective country. Minimum requirements: 3 min. at 134°C (EU: 5 min. at 134°C) Maximum sterilization temperature: 138° C

- Drying time: For steam sterilization, we recommend a drying time of 15 to 40 minutes. Choose a suitable drying time depending on the autoclave and load. Refer to the autoclave's instructions for use. The use of steam sterilization, strong cleaning and disinfection fluids (alkaline pH>9 or acidic pH<5) can reduce the service life of appliances. The manufacturer declines all responsibility in these cases. Appliances should not be exposed to temperatures above 136°C. The user has an obligation to ensure that the reprocessing processes, including resources, materials and personnel, are capable of achieving the required results. State-of-the-art and often national legislation require that these processes and the resources included be properly validated and maintained.

Note: We reserve the right to modify the equipment design, technology, accessories, instruction manual, and contents of the original packing list at any time without notice. If there are differences between the model and the actual equipment, the actual equipment should be considered the norm.

After sterilization:

2. Remove the product from the autoclave.
6. Allow the product to cool to room temperature for at least 30 minutes. Do not use additional cooling. Check that sterilization wrappers or bags are not damaged.

Rapid sterilization is not allowed for products with lumen.

The manufacturer is not responsible for the use of other sterilization processes (e.g., ethylene oxide, formaldehyde, and low-temperature plasma sterilization). In these cases, it is necessary to observe the respective valid standards (EN ISO 14937/ANSI AAMI ISO 14937 or the procedure-specific standard) and check the suitability and efficacy in principle of the procedure (if necessary, including investigations on residues of sterilizing agents), taking into account the specific geometry of the product as part of the validation. Maximum sterilization temperature 138°C

15. Storage: Store sterilized products in a dry, clean and dust-free environment, with a relative humidity of 10% to 93%, an atmospheric pressure of 70kPa to 106kPa and a temperature of -20°C to +55°C; Refer to the label and instructions for use. After sterilization, the product should be packed in a medical sterilization bag or clean resealable container and stored in a special storage cabinet. The storage time should not exceed 7 days. If it is exceeded, it must be reprocessed before use.

16. Life shelf: The products are designed for a large number of sterilization cycles. The materials used in its manufacture have been selected accordingly. However, with each new preparation for use, thermal and chemical stresses will cause the devices to age. If the number of permitted re-sterilization cycles is limited, this will be indicated in the product-specific instructions.

Shelf life: 5 years

The use of ultrasonic baths and strong cleaning and disinfection liquids (alkaline pH>9 or acidic pH<5) may reduce the service life of the devices. The manufacturer declines all responsibility in these cases.

The devices must not be exposed to temperatures above 138°C. It is the user's obligation to ensure that the reprocessing processes, including resources, materials and personnel, are capable of achieving the required results. State of the art and often national legislation requires that these processes and the resources included are properly validated and maintained.

Note: We reserve the right to modify the equipment design, technology, accessories, operating instructions and contents of the original packing list at any time and without prior notice. If there are differences between the model and the actual equipment, the actual equipment should be considered the standard.

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	Caution
	Refer to instruction manual/booklet
	Manufacturer
	Date of manufacture
	Handle with care
	Keep dry
	For indoor use only
	Can be autoclaved
	Recovery
	Temperature limitation -20° C-+40° C
	Humidity limitation 10%-93%
	lot number
	Medical device
	CE marking with identification number of the Notified Body
	Authorized representative in the European Community
	Green: No Wear-Tip is OK Tip replacement is not necessary
	Yellow: Wear of 1mm-Tip is showing some wear Tip replacement is not necessary
	Red: Wear of 2 mm-Tip is badly worn Tip replacement is necessary
	Effective length

				
ECO DESIGN	CLOSE MANUFACTURE	RESPONSIBLE CONSUMPTION	SMART WASTE	ONLINE DISTRIBUTION

	
Guilin Refine Medical Instrument Co. Ltd. No.8-3, Information IndustrialPark,High Tech Zone,Qixing District, Guilin, Guangxi 541004, P.R.China	Dentared Odontology Services SL C. Torres Quevedo 13 Alcobendas, 28108, Madrid, Spain
	
Med Net EC-Rep GmbH Borkstrasse 10, 48163 Münster, Germany	www.vidu.eco

					
CE	0123	MD	No sharp	134°C	Information