Vinyl polysiloxane impression material  
DIN EN ISO 4823



# PANASIL<sup>®</sup>

**INSTRUCTIONS FOR USE** / Gebrauchsinformation / Mode d'emploi / Gebruiksaanwijzing /  
Инструкция по использованию / Istruzioni per l'uso / Instrucciones de uso / Kullanım Kılavuzu /  
Instruções de uso / Instrukcja użycia / Інструкція для застосування / 取扱説明書 / تعليمات الاستخدام

# CONTENT

Inhalt / Contenu / Inhoud / содержимое /  
Contenuto / Contenido / İçerik / Conteúdo /  
Treść / Зміст / 目次 / الفهرس

|                  |      |
|------------------|------|
| Deutsch .....    | 4    |
| English .....    | 16   |
| Français .....   | 27   |
| Nederlands ..... | 39   |
| Русский .....    | 51   |
| Italiano .....   | 64   |
| Español .....    | 76   |
| Türkçe .....     | 88   |
| Português .....  | 100  |
| Polski .....     | 112  |
| Українська ..... | 126  |
| 日本語 .....        | 140  |
| 153 .....        | عربي |

Vinyl polysiloxane impression material  
DIN EN ISO 4823

**Panasil®** Putty  
**Panasil®** Putty Fast  
**Panasil®** Putty Soft  
**Panasil®** binetics Putty Fast  
**Panasil®** binetics Putty Soft  
**Panasil®** tray Fast Heavy  
**Panasil®** tray Soft Heavy  
**Panasil®** tray Soft Heavy Fast

**Panasil®** monophase Medium  
**Panasil®** initial contact Regular  
**Panasil®** initial contact Light  
**Panasil®** initial contact Light Fast  
**Panasil®** initial contact X-Light  
**Panasil®** initial contact X-Light Fast  
**Panasil®** contact two in one Light  
**Panasil®** contact plus X-Light

**For dental professional use only.**

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## Product Description

Panasil® impression materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation. Panasil® initial contact displays a high initial hydrophilicity.

The Panasil® family consists of five different viscosities (putty, heavy bodied, medium bodied, light bodied, x-light bodied), available in an assortment of delivery systems: standard 1:1 (50 ml automix cartridges) and 5:1 (380 ml cartridges) for use in most commercially available other automatic dispensing and mixing machines (e.g., Sympress), and traditional 1:1 putty jars.

## Indications / Techniques

Panasil® Putty (Putty, Putty Fast, Putty Soft) and Panasil® binetics Putty (Putty Fast, Putty Soft) are to be used as preliminary materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step putty-wash impression technique using a foil (plastic putty spacer)

- One-step putty impression technique for forming functional peripheries

Panasil® tray Fast Heavy, Panasil® tray Soft Heavy and Panasil® tray Soft Heavy Fast are to be used as heavy bodied materials for:

- One-step impression technique (simultaneous technique) using single or dual viscosities
- Two-step impression technique using dual viscosities (Panasil® tray Fast Heavy)
- Functional impressions (one-step: Panasil® tray Soft Heavy)

Panasil® monophasic Medium is to be used as a medium bodied tray or syringeable impression material for:

- Taking impressions over fixed/removable restorations and implants (i.e., transferring impression posts and bridge components)
- Functional impressions
- Fabricating crown and bridgework or inlays
- Fabricating full or partial dentures
- Reline impressions

- Use in the simultaneous mixing technique as well as the putty-wash and Multi tray techniques
- Transferring root posts when fabricating posts and cores indirectly

Panasil® initial contact, Panasil® contact plus X-Light and Panasil® contact two in one Light are to be used as syringeable impression materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step impression technique using a spacer
- One-step impression technique (simultaneous technique) using dual viscosities
- Reline impressions
- Fabricating full or partial dentures

## Patient Target Group

Patients being treated in the context of a dental procedure.

## Warnings

Do not use Panasil® impression materials as a temporary reliner.

Do not use Panasil® impression materials with condensation-curing silicones, polyether, Vinylsiloxanether® or polysulfide materials.

Panasil® Putty impression materials are not suitable for detailed impressions when used alone.

## Potential Allergenic Substances

Platinum compound, 4-methoxyphenol.

## Cautions

Do not use after expiration date.

Normally, impression materials do not cause allergic reactions; however, an allergic reaction cannot be ruled out in sensitive individuals.

If in doubt, we recommend performing an allergy test before use.

Do not use the material in case of known allergies to any of the ingredients or in case of contact allergies.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate,

polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

Do not interchange the base paste and catalyst paste lids or scoops for hand-mixed jar putty.

To ensure an optimal connection between the two materials, after the trial positioning, the preliminary impression should also be rinsed with water and dried before the correction material is applied. Both steps of the putty-wash impression technique (preliminary impression and wash impression) should also be carried out in quick succession.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth/jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

## Note

To ensure optimal impressions, the temperature of the material should not deviate from 23°C (74°F) before applying. Otherwise, working and setting times will be affected.

## Impression tray: Preparation and Adhesives

In general, all conventional impression trays can be used if the necessary pressure built up can be generated. If it is not possible to achieve sufficient retention with the impression material, brush the impression tray with Panasil® adhesive. The surface of the impression tray must be repeatedly and fully brushed so that the entire surface of the tray is visibly blue. Allow the adhesive to dry in accordance with the information in the direction for use and cover with sufficient impression material in the impression tray.

## Instructions for Use: Jar Material for Hand-mix Putty Preparation

Only use materials with the same lot number.

Do not interchange the base paste and catalyst paste lids or scoops for hand-mixed jar putty.

Use the scoops to dispense equal amounts of base paste and catalyst paste.

The different color scoops must only be used for the materials with corresponding colors.

Close the containers carefully after use, and ensure that the lids are not interchanged.

Contamination of base paste and catalyst paste in the container renders the material unusable.

Knead the base paste and catalyst paste for 45 seconds until the material is homogeneously colored.

Should gloves be worn, test them for compatibility prior to mixing with a sample of the material to be mixed. Certain types of gloves can prevent the polymerization (e.g., latex gloves) or impair the adhesive bond (e.g., nitrile gloves). Using vinyl gloves is recommended.

Prior to taking the wash impression, the preliminary impression must be cleaned with water, dried and trimmed using standard methods.

## **Instructions for the Safe Handling of the Kettenbach 380-mL cartridge system; 5:1**

### **Components: 380-mL cartridge, Dynamic Mixers**

Prior to first use please unscrew the blue plug counterclockwise and pull it out (figure 1, page 164). Check the outlet opening. Then screw on the dynamic mixer (figure 2, page 164). Ensure that the positioning aid is on the top.

Please use only the dispensing and mixing components compatible with the Kettenbach 380 ml cartridge system for dispensing and mixing impression material.

The Kettenbach 380 ml cartridge is heavy and could be damaged if dropped. Damaged cartridges must not be used under any circumstances.

When using the material in a Sympress II dispensing device, please activate the speed setting indicated on the product label (figure 6, page 165). Failure to comply with this instruction may modify the curing characteristics of the impression materials.

If you are using a unit other than the Sympress mixing device, please follow the instructions for the automatic dispensing and mixing system you are using.

When changing to a different mixing device, always discard at least 3 cm of mixed material.

## **Using the Kettenbach 380-mL cartridge system; 5:1 in Dispensing Device Sympress / Sympress II**

By pressing the plunger return knob, allow the plunger to return to its original position automatically (while the cover is closed; figure 3, page 164). Next, open the cover (figure 4, page 164).

Insert the Kettenbach 380 ml cartridge with the pre-assembled mixer in the Sympress/Sympress II dispensing and mixing unit (figure 5, page 165). If you are working with a previously used cartridge, the filled mixer will already be located on the cartridge as a cap.

If you are using a Sympress II dispensing device, please activate the speed setting indicated on the product label. A light signal will confirm your selection (figure 6, page 165).

Close the device cover (figure 7, page 165).

## Using the Cartridge for the First Time

Dispense the impression material (ensure that the dynamic mixer is wound into place fully) until the base material and catalyst are mixed homogeneously; discard at least the first 3 cm when dispensing (figure 8, page 165).

At least 3 cm of the mixed material must be discarded each time the plunger is restarted or when using the cartridge in different mixing units.

Load the impression tray (or application syringe) with the required amount of material (figure 9, page 165) by pressing the appropriate feed key. Hold the impression at an angle when loading and press against it lightly. Leave the mixer in the material when dispensing to avoid air bubbles.

Release the feed key as soon as the required amount of material has been dispensed. The filled dynamic mixer should remain on the cartridge as a seal.

Remove the used mixer before the next application by turning it to the left (figure 10, page 165), visually check the outlet opening for blockage and attach a new dynamic mixer (figure 10, page 165).

## Further Use of a Cartridge Already in Use

Wind a new dynamic mixer onto the same mixing unit and begin the new application (material no longer has to be discarded).

## Instructions on the Use of Cartridge Material in the Applyfix 4 Dispensing Gun DS-50 1:1/2:1

Insert the cartridge into the Applyfix 4 dispensing gun DS-50 1:1/2:1.

Ensure that the notches on the cartridge base are pointing downwards. The clasp will not close if the cartridge is not properly inserted into the dispensing gun.

The cartridge's tamper-evident closure can now be turned (see Fig. 1).


Abb. 1: Opening the tamper-evident closure

The tamper-evident closure ensures that this is an intact original product from Kettenbach Dental.

When used for the first time, the cartridge closure is visibly and audibly separated from the cartridge body. After use, the closure cap can be used to close the cartridge.

Insert the mixing tip with the integrated positioning aid into the coding hole provided for that purpose and rotate clockwise into the thread (see Fig. 2).


Abb. 2: Inserting the mixing tip

The intuitive handling significantly reduces mistakes when inserting.

Push the plunger into the cartridge and begin applying immediately.

The mixing system ensures even delivery at all times, so it is no longer necessary to throw out the impression material.

Either dispense the necessary amount of impression material directly into the impression tray, on to the preliminary impression or into an impression syringe. Use the cartridge with the mixer attached and intraoral tip or the impression syringe to apply the material around the prepared teeth.

Replace the closure cap after use. Before using the cartridge again remove the closure cap, attach a new mixing tip and proceed as described.

## Cleaning and Disinfection of the Applyfix 4

Remove the cartridge.

Remove any uncured paste using a cloth soaked in alcohol.

The Applyfix 4 dispensing gun can be used several times if properly disinfected. Kettenbach recommends exchanging the gun if there are visible traces of use. Disinfection can be carried out using a commercially available dipping disinfectant. Kettenbach recommends the use of glutaraldehyde-based solutions.

## Storage Conditions

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

## Disinfection

Impressions can be disinfected using, for example, a 2 % acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials. See manufacturer instructions.

## Model Preparation

Impressions with Panasil® may be poured immediately after disinfection of the impression up to several weeks with high quality dental stone (class III–V). Impressions can also be poured several times.

## Composition

Vinyl polysiloxanes, Hydrogen polysiloxanes, Filler, Pigments, Additives, Platinum compound.

## Disposal

The empty cartridges / mixing tips / applicator tips and set material can be disposed in the household rubbish.

## Trademarks

Panasil® is a registered trademark of Kettenbach GmbH & Co. KG.

## Note

Selected Kettenbach materials are available in certain markets only.

For dental professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

If serious incidents occur that are related to the product, please immediately contact the manufacturer and the competent authority.

## Limitation of Liability

Except where prohibited by law, Kettenbach GmbH & Co. KG will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

The information provided for Kettenbach products is based on comprehensive research and experience in application technology. Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development. However, users must comply with and consider all recommendations and information in connection with any use.

Updated: May 10, 2023

## Product Overview

| Product name                  | ISO 4823             | Consistency (approx.) mm | Mixing ratio and total content volume per unit | Mixing technique                                                        | Mixing element             |
|-------------------------------|----------------------|--------------------------|------------------------------------------------|-------------------------------------------------------------------------|----------------------------|
| Panasil® Putty                | Type 0, Putty        | 24                       | 1:1, 900 ml in jars                            | Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste) | Hand mix                   |
| Panasil® Putty Fast           | Type 0, Putty        | 24                       | 1:1, 900 ml in jars                            | Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste) | Hand mix                   |
| Panasil® Putty Soft           | Type 0, Putty        | 24                       | 1:1, 900 ml in jars                            | Manual, 1 scoop (12.3 ml) per component (base paste and catalyst paste) | Hand mix                   |
| Panasil® binetics Putty Fast  | Type 0, Putty        | 23                       | 5:1, 380 ml cartridge                          | Sympress/Sympress II or other automatic dispensing and mixing unit      | Dynamic mixer, blue        |
| Panasil® binetics Putty Soft  | Type 0, Putty        | 23                       | 5:1, 380 ml cartridge                          | Sympress/Sympress II or other automatic dispensing and mixing unit      | Dynamic mixer, blue        |
| Panasil® tray Fast Heavy      | Type 1, Heavy bodied | 32                       | 5:1, 380 ml cartridge                          | Sympress/Sympress II or other automatic dispensing and mixing unit      | Dynamic mixer, blue        |
| Panasil® tray Fast Heavy      | Type 1, Heavy bodied | 34                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                               | Mixing tip, blue, Ø 6.0 mm |
| Panasil® tray Soft Heavy      | Type 1, Heavy bodied | 32                       | 5:1, 380 ml cartridge                          | Sympress/Sympress II or other automatic dispensing and mixing unit      | Dynamic mixer, blue        |
| Panasil® tray Soft Heavy      | Type 1, Heavy bodied | 34                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                               | Mixing tip, blue, Ø 6.0 mm |
| Panasil® tray Soft Heavy Fast | Type 1, Heavy bodied | 34                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                               | Mixing tip, blue, Ø 6.0 mm |

## Product Overview

| Product name                          | ISO 4823              | Consistency (approx.) mm | Mixing ratio and total content volume per unit | Mixing technique                                                     | Mixing element             |
|---------------------------------------|-----------------------|--------------------------|------------------------------------------------|----------------------------------------------------------------------|----------------------------|
| Panasil® monophase Medium             | Type 2, Medium bodied | 35                       | 5:1, 380 ml cartridge                          | Sympress / Sympress II or other automatic dispensing and mixing unit | Dynamic mixer, blue        |
| Panasil® monophase Medium             | Type 2, Medium bodied | 35                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, blue, Ø 6.0 mm |
| Panasil® initial contact Regular      | Type 2, Medium bodied | 37                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, blue, Ø 6.0 mm |
| Panasil® initial contact Light        | Type 3, Light bodied  | 41                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, red, Ø 4.0 mm  |
| Panasil® initial contact Light Fast   | Type 3, Light bodied  | 41                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, red, Ø 4.0 mm  |
| Panasil® initial contact X-Light      | Type 3, Light bodied  | 44                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, red, Ø 4.0 mm  |
| Panasil® initial contact X-Light Fast | Type 3, Light bodied  | 44                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, red, Ø 4.0 mm  |
| Panasil® contact two in one Light     | Type 3, Light bodied  | 38                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, blue, Ø 6.0 mm |
| Panasil® contact plus X-Light         | Type 3, Light bodied  | 42                       | 1:1, 50 ml cartridge                           | Applyfix 4 dispensing gun DS-50 1:1 / 2:1                            | Mixing tip, red, Ø 4.0 mm  |

## Technical Data

| Product name                  | Mixing ratio | Working time at 23°C/ 74°F ≤ | Working time at 35°C/ 95°F ≤ | Intraoral setting time at 35°C/ 95°F ≥ | Total setting time* ≥ | Hardness (approx.) Shore | Linear dimensional change (approx.) % | Elastic Recovery test (approx.) % | Strain in compression (approx.) % |
|-------------------------------|--------------|------------------------------|------------------------------|----------------------------------------|-----------------------|--------------------------|---------------------------------------|-----------------------------------|-----------------------------------|
| Panasil® Putty                | 1:1          | 2 minutes                    | not applicable               | 2 minutes                              | 4 minutes             | A 66                     | - 0.20                                | 99.0                              | 2.7                               |
| Panasil® Putty Fast           | 1:1          | 1 minute 30 seconds          | not applicable               | 2 minutes                              | 3 minutes 30 seconds  | A 66                     | - 0.20                                | 99.0                              | 2.7                               |
| Panasil® Putty Soft           | 1:1          | 2 minutes                    | not applicable               | 2 minutes                              | 4 minutes             | A 60                     | - 0.20                                | 99.0                              | 2.7                               |
| Panasil® binetics Putty Fast  | 5:1          | 1 minute 30 seconds          | not applicable               | 2 minutes 30 seconds                   | 4 minutes             | A 63                     | - 0.20                                | 99.5                              | 2.5                               |
| Panasil® binetics Putty Soft  | 5:1          | 2 minutes                    | not applicable               | 3 minutes                              | 5 minutes             | A 56                     | - 0.20                                | 99.5                              | 3.5                               |
| Panasil® tray Fast Heavy      | 5:1          | 1 minute 20 seconds          | not applicable               | 2 minutes                              | 3 minutes 20 seconds  | A 62                     | - 0.20                                | 99.7                              | 2.5                               |
| Panasil® tray Fast Heavy      | 1:1          | 1 minute                     | not applicable               | 2 minutes                              | 3 minutes             | A 62                     | - 0.20                                | 99.5                              | 3.0                               |
| Panasil® tray Soft Heavy      | 5:1          | 2 minutes                    | not applicable               | 2 minutes                              | 4 minutes             | A 55                     | - 0.20                                | 99.7                              | 3.0                               |
| Panasil® tray Soft Heavy      | 1:1          | 2 minutes                    | not applicable               | 2 minutes                              | 4 minutes             | A 55                     | - 0.20                                | 99.5                              | 3.0                               |
| Panasil® tray Soft Heavy Fast | 1:1          | 1 minute                     | not applicable               | 2 minutes                              | 3 minutes             | A 55                     | - 0.20                                | 99.5                              | 3.0                               |

## Technical Data

| Product name                          | Mixing ratio | Working time at 23°C/ 74°F ≤ | Working time at 35°C/ 95°F ≤ | Intraoral setting time at 35°C/ 95°F ≥ | Total setting time* ≥ | Hardness (approx.) Shore | Linear dimensional change (approx.) % | Elastic Recovery test (approx.) % | Strain in compression (approx.) % |
|---------------------------------------|--------------|------------------------------|------------------------------|----------------------------------------|-----------------------|--------------------------|---------------------------------------|-----------------------------------|-----------------------------------|
| Panasil® monophase Medium             | 5:1          | 2 minutes                    | 1 minute                     | 2 minutes                              | 4 minutes             | A 60                     | – 0.20                                | 99.7                              | 3.5                               |
| Panasil® monophase Medium             | 1:1          | 2 minutes                    | 1 minute                     | 2 minutes                              | 4 minutes             | A 60                     | – 0.20                                | 99.7                              | 3.5                               |
| Panasil® initial contact Regular      | 1:1          | 1 minute 30 seconds          | 1 minute                     | 2 minutes 30 seconds                   | 4 minutes             | A 46                     | – 0.20                                | 99.7                              | 3.0                               |
| Panasil® initial contact Light        | 1:1          | 1 minute 30 seconds          | 1 minute                     | 2 minutes 30 seconds                   | 4 minutes             | A 46                     | – 0.20                                | 99.3                              | 3.5                               |
| Panasil® initial contact Light Fast   | 1:1          | 1 minute                     | 30 seconds                   | 2 minutes                              | 3 minutes             | A 46                     | – 0.20                                | 99.3                              | 3.5                               |
| Panasil® initial contact X-Light      | 1:1          | 1 minute 30 seconds          | 1 minute                     | 2 minutes 30 seconds                   | 4 minutes             | A 46                     | – 0.20                                | 99.3                              | 3.5                               |
| Panasil® initial contact X-Light Fast | 1:1          | 1 minute                     | 30 seconds                   | 2 minutes                              | 3 minutes             | A 46                     | – 0.20                                | 99.3                              | 3.5                               |
| Panasil® contact two in one Light     | 1:1          | 2 minutes                    | 1 minute                     | 2 minutes                              |                       | A 46                     | – 0.20                                | 99.7                              | 3.5                               |
| Panasil® contact plus X-Light         | 1:1          | 2 minutes                    | 1 minute                     | 2 minutes                              |                       | A 46                     | – 0.20                                | 99.7                              | 5.0                               |

\* Total setting time (removal time from mouth) from start of mix

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