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1. DESCRIPCIÓN

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INSTRUCTIONS FOR USE:

1. DESCRIPTION

VENTURA NG CAP belongs to the latest generation of ternary alloys with a high copper content and a mixture of spherical particles; being the latter predominantly present. Its balanced formula and ultra-fine granulometry enable amalgams to be prepared without the presence of the gamma 2 phase and with very high compression resistance. **VENTURA NG CAP** predosed self-activating capsules contain **VENTURA NG** alloy in powder and tridistilled mercury in the proportion indicated in the table.

2. INDICATIONS

Fillings that must stand pressure in the back teeth (Class I and II).

3. CONTRAINDICATIONS

The use of amalgam is not recommended when the doctor considers it unsuitable for the patient. It is not recommended to make or extract a filling of dental amalgam in patients who are hyper-sensitive to mercury. Avoid using dental amalgam in children under fifteen years old. It is not recommended to make or extract a filling of dental amalgam in women who are pregnant or breast-feeding. Avoid contact between dental amalgam and other existing metallic restorations.

4. MERCURY - ALLOY PROPORTION

VENTURA NG CAP predosed capsules contain mercury and alloy in a 0,94:1 proportion.

5. AMALGAMATION

VENTURA NG CAP capsules must be processed 7 ± 1 seconds on a high energy amalgamator (>4200 rpm) such as ventura mix 2.

	Capsule	Alloy	Hg
VENTURA NG CAP 1	RED SKY BLUE	400 mg	377 mg
VENTURA NG CAP 2	RED BLUE	600 mg	565 mg
VENTURA NG CAP 3	RED WHITE	800 mg	753 mg

6. PROPERTIES

VENTURA NG CAP fully fulfil the requirements set by standards ISO 24234, ISO 20749 and ISO 23325. Once the capsules have undergone vibration for the ideal length of time a shiny ball of amalgam is produced, which comes cleanly away from the capsule.

7. INFORMATION OF INTEREST

Drilling, polishing and grinding of fillings must always be accompanied by cooling with water and vacuum suction. Mercury reacts with and embrittles particular metals and their alloys. Avoid unnecessary contact between mercury and those metals (and their alloys). If moisture is introduced into the amalgam before it has set, properties such as strength and corrosion resistance may be affected adversely. If the alloy contains zinc, such contamination may result in an excessive expansion (delayed expansion). Whenever it is possible, use a dry field.

8. RECOMMENDATIONS FOR SAFETY PRECAUTIONS

Mercury presents a health hazard if handled incorrectly. Mercury is toxic by vapor inhalation and the effect is cumulative. Spillages of mercury should be removed immediately, including from places which are difficult to access. Use a plastic syringe to draw it up. Smaller quantities can be covered by sulfur powder and removed. Individual small droplets can be picked up by tin foil (Sn) and removed. During and after cleanup thoroughly ventilate the area where the spill has occurred. Vacuum cleaners should not be used. Avoid inhalation of mercury vapours. Keep container tightly closed. Keep locked up and out of reach of children. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Very toxic to aquatic organisms. It may cause long-term adverse effects in the aquatic environment. Waste material and all primary containers that have held mercury shall be disposed of following appropriate management practice.

This information is in compliance with ISO 24234 and ISO 20749.

9. RESTRICTIONS TO BE APPLIED WITHIN THE EU ACCORDING TO REGULATION 2017/852

From 1 July 2018, dental amalgam shall not be used:

- for dental treatment of deciduous teeth,
- for children under 15 years,
- for pregnant or breastfeeding women, except when deemed strictly necessary by the dental practitioner based on the specific medical needs of the patient.

Material developed for the exclusive use by dental professionals.

Instructions for use issued in October 2021 (2021-10).

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