R-020195

ESINYECTORA DE RESINAS	2
ENRESIN INJECTION SYSTEM	5
FRSYSTEME D'INJECTION DE RESINE	9

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

Rev. 28/06/04

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This system of introduction of the resin in the flask for the manufacture of prosthesis by injection is an operating and practical system, without forcing the operator to change:

- The brands and products you usually use.
- The usual technique of preparation of the resin.
- Technology and devices for polymerization.

It avoids the acquisition of particularly complex and expensive equipment.

SYSTEM COMPONENTS

1. Flsk for injection: made of hot stamped and chromed brass.
2. Flange for injection, in 3 versions, for:
 - 1 flask.
 - 2 flasks.
 - 3 flasks.

built in forged brass, very robust, work without difficulty, with or without resin inside.

3. Syringe for injection: built in stainless steel, for easy cleaning and maintenance. It is equipped with a key to screw the syringe and the stopper of the flask. It is supplied with replacement piston rubber ring, which needs periodic replacement (can be purchased in oleopneumatic establishments).

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INSTRUCTIONS FOR USE

In the case of total prosthesis, it is advisable to position the prosthesis by placing the retromolar tubers (posterior part of the last molar) in the part of the threaded hole to insert the syringe and the anterior part of the prosthesis positioned on the opposite side, where, in the perimeter of the flask, there are two holes of 2 mm in diameter, which serve for the exit of the excess resin.

1. After the melting of the denture wax in a flask, as is normally the case, 1 or 2 channels of approximately 2 mm are to be excavated in the plaster, connecting the anterior gingival part with the injection discharge holes, which are exactly in the same part of the flask (in effect we have equipped the flask with 4 discharge holes; with time and experience, excellent results will be obtained, both in the positioning of the prosthesis and in the increase in the number of casts, achieving excellent results in prostheses particularly large or complicated).

A 2 mm cutter or cutter can be used for this operation.

2. Now dig in the posterior part of the prosthesis another 2 channels of diameter superior to the previous ones that connect the retromolar tubers with the injection hole, both converging in the single main hole of the flask, which will serve to lodge the syringe.
3. Clean plaster waste thoroughly. Then apply the usual insulation.
4. Close the two half-fillets and clamp in the longitudinal direction, so that the two flange columns meet at the apex of the flask and do not obstruct either the inlet hole or the resin outlet.
5. Tighten the flange decisively, but without exaggerating, not hammering on the tightening handle. This way of working eliminates the possibility of some prominence or elevation of chewing.
6. Unscrew the front cap of the syringe (needle side) and screw it into the flask through the narrowest part. Push the piston back into the syringe barrel until it stops, unscrewing the threaded rod by means of the transverse handle counterclockwise. In this way the body of the syringe remains completely empty and free to be filled with resin. Note that the capacity of the syringe allows filling at least 2 flasks of average size, one lower and one higher.

ATTENTION: the plasticity of the resin should be slightly softer than that normally used

7. Once the syringe has been filled with resin, it is screwed into the cap previously threaded in the flask.
8. Now begin to introduce the resin by threading the handle of the syringe in a clockwise direction. Perform the operation PROGRESSIVELY, without return, so as not to produce pores or air bubbles in the resin. The injection phase is completed when two compact resin threads come out through the 2 discharge holes.

9. At this moment, the pressure of the resin inside the flask is such that it not only completely fills any cavity, but also makes it difficult to unscrew the syringe from the flask; for this reason it has been equipped with a key that allows, without wasting time, unscrewing the syringe and filling the remaining flanged flasks.
10. When the syringe is removed from the thread, the resin will tend to come out through the same hole. Now we must insert the stopper of the flask in the place of the syringe. While screwing thoroughly you will see that the resin threads continue to come out from the opposite side of the flask, indicating this a perfect execution of the filling operation.
11. Proceed to polymerize and execute all successive operations as usual. Once the finished work, you can see a particularly clean prosthesis, very easy to deburr, polish and polish, without forgetting the undoubted improvements obtained in the physical characteristics of the prosthesis due to molecular compaction and the absence of porosities or savings substantial time. In addition, we want to emphasize the absolute elimination of the annoying prominences or elevations of mastication

OPERATING PRINCIPLE

The pressure exerted by the piston inside the syringe pushes the resin gradually into ever narrower sections, so that any cavity in the hard plaster is necessarily filled, before the resin comes to discharge the pressure through the holes of download.

In short, there is a continuous and gradual decrease in the section through which the resin passes, according to the following measures:

Inner diameter of the syringe barrel.....	24 mm
Diameter in the cap of the siringe	4 mm
Diameter channels from syringe to prosthesis	2 x 2.5 mm
Diameter channels discharge behind the prosthesis.....	2 x 2.5 mm

TIPS

- If a plaster / plaster mixture is normally used, it is not necessary to change the habits, although it should be taken into account that using gypsum lasts an exceptional finish.
- Use any type of acrylic or methacrylic resin to polymerize hot, preferably the one that is known and handled well.
- Number of casting and unloading channels to be made:
 - In partial or total inferior prostheses: 1 injection and 1 discharge, making the entry on a retromolar and the exit at the opposite end of the prosthesis.
 - In complete lower prostheses with a particularly high gingival crest, use a traditional system, since in some area a bubble may form in which the air is not replaced by the resin. Only with the experience that gives the practice can these cases be overcome, on the other hand, rare if soft resin is used.

NOTE: In any case, this system is an integral evolution of the traditional system and therefore allows to fill the flasks of resin by the classical method when you doubt the perfect filling of the flask (opening the flask, removing the syringe and placing the cap in its place).

MAINTENANCE

The flange does not require special maintenance:

English

- ◆ In the flask it is sufficient to keep the threaded inlet hole and discharge holes clean.
- ◆ Very important is the cleaning of the syringe that must always be done before the resin hardens too much and makes that operation difficult.
- ◆ Completely disassemble the syringe and carefully wash all the parts with hot water; particularly remove the black rubber ring from its housing in the piston and remove any resin residue from its accommodation throat.
- ◆ Lightly lubricate the threaded rod of the piston so that it slides better.

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