· Conscious dentistry ·

VIDU


EN

ADDITION CURING SILICONE Light Body


Polyvinylsiloxane for high precision dental impressions, addition-curing. Product complies with ISO 4823. For professional use. Rx ONLY.

INSTRUCTIONS FOR USE

APPLICATIONS: Silicone indicated for the two-step putty/wash and for the sandwich technique. The impression can be electroplated.

PREPARATION: Dry the areas to be reproduced and prevent further inflow of liquid into them. Load the cartridge into the dispenser (1) and after removing the cap, extrude a small amount of silicone from the two nozzles, moving the lever several times (2). Insert the mixing tip into the cartridge and if necessary, place the small oral tip on this for perfect and precise gingival extrusion. Actuate the lever forcefully so as to extrude the ADDITION CURING SILICONE Light Body VIDU from the mixing tip, loading the mixed material into the impression tray (3) and insert it into the mouth (4). It is important to check that the extruded silicone is homogeneously coloured. To preserve any product remaining in the cartridge correctly, leave the mixing tip fixed until the next use. Before applying a new mixing tip, ensure that the silicone comes out of the nozzles regularly. If the cartridge remains unused for a long period, it may become blocked by hardened silicone, which must be removed with a thin instrument before a new mixing tip is placed on.

TECHNICAL DATA	SUPERLIGHT		LIGHT	
	NORMAL	FAST	NORMAL	FAST
ISO 4823 classification	Type 3		Type 3	
Mixing time	(automixing)		(automixing)	
Working time	2'	120"	2'	120"
Minimum time in mouth	130"	60"	130"	60"
Recovery from deformation	99,9%		99,9%	
Strain in compression	4%		4%	
Linear dimensional change at 24 hours	- 0,02%		- 0,02%	
Detail reproduction	5µm		5µm	
Hardness Shore A	47		47	

*Measurement at 23 +/- 2°C and 50 +/- 10% relative humidity.

WARNINGS: The above times can be influenced by:
 - variations in the temperature of the material and/or environment (higher temperatures reduce them while lower temperatures prolong them);
 - humidity.

Avoid direct or indirect contact with latex rubber dams, gloves or with components of such material. Contact with latex inhibits the hardening of polyvinylsiloxanes. Accurately wash your hands and use vinyl/nitrile gloves.

To eliminate traces of saliva or blood after taking the impression, rinse under running water and dry well. The impression can be stored without any special precautions. Cast the model in plaster at least 30 minutes after removal from the mouth.

For disinfection, immerse the impression for 5 minutes in a solution of glutaraldehyde (concentration 0.5%-2%) or sodium hypochlorite (concentration 0.1%-1%). Rinse the impression under running water after disinfection.

IMPRESSION TRAY: All types of rigid impression trays can be used.
GUARANTEE: The manufacturer guarantees the quality of its products and undertakes to replace them whenever defects in production are found within the expiry dates. The manufacturer will not be held responsible and will not recognize any damage that is direct, indirect, and accidental or caused by incorrect use of the material. Check the suitability of the product prior to use.

STORAGE: Shelf life 3 years when stored a 77 °F (25 °C).
PRESENTATION: - [COD V152013] ADDITION CURING SILICONE Light Body, Normal Set.
 - [COD V152014] ADDITION CURING SILICONE Light Body, Fast Set.

ADDITION CURING SILICONE Light Body

ES · SILICONA DE CURADO POR ADICIÓN Light Body
IT · SILICONE A POLIMERIZZAZIONE PER ADDIZIONE Light Body
FR · SILICONE PAR ADDITION A POLYMERISATION Light Body
PT · SILICONE POLIMERIZÁVEL POR ADIÇÃO Light Body
DE · ADDITIONSVERNETHZENDES SILIKON Light Body

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