ADVANCED ENDODONTIC  
WHITENING SYSTEM  
INSTRUCTIONS FOR USE

Pola Professional 37.5% is a hydrogen peroxide whitening system designed to medically bleach endodontically treated teeth.

**INDICATIONS:**

- Whitening of discoloured non-vital teeth that have been endodontically treated.

**COMPOSITION:** Contains 37.5% Hydrogen peroxide.

**NET WEIGHT:**

Depending on kit purchased. Please refer to outside of box:

- 3.3g Pola Professional 37.5% syringe
- 1g Gingival Barrier

**SUGGESTIONS FOR SAFER NON-VITAL BLEACHING:**

**1. Verify adequate endodontic obturation.**

The root canal filling must have been completed to appropriate endodontic standards. The root canal should be completed at least a month before bleaching. The quality of root canal obturation should always be assessed clinically and radiographically prior to bleaching. Adequate obturation ensures a better overall prognosis of the treated tooth. It provides an additional barrier against damage by peroxide to the periodontal ligament and periapical tissues.

**2. Use protective barrier.**

This is essential to prevent leakage of bleaching agents which may infiltrate between the gutta-percha and root canal walls, reaching the periodontal ligament via dentinal tubules, lateral canals or the root apex. The ideal barrier should protect the dentinal tubules.

**Caution:** Seal gutta-percha to a level which is coronal to potential communication with periodontal tissue. We recommend using a glass ionomer cement such as Riva SC (SDI Ltd).

**3. Avoid acidic etching.**

Acid etching has not proven to be beneficial (Casey et al 1989) and the use of caustic chemicals in the pulp chamber is undesirable, as periodontal ligament irritation may result.

**4. Recall patients periodically.**

Bleached teeth should be frequently examined both clinically and radiographically. Root resorption may occasionally be detected as early as 6 months after bleaching. Early detection improves the prognosis as corrective therapy may still be applied.

**TIPS:**

- Isolate tooth effectively. Use of a rubber dam is advised.
- Protect oral mucosa to prevent chemical burns.
- Advise the patient that a whitened non-vital tooth will not look the same as a vital tooth.
- Delay any esthetic restorations for 2 to 4 weeks after treatment.
- Advise patient not to smoke, consume highly coloured, hot or cold foods or beverages for at least 2 days.
- Recall patients periodically to examine the treated teeth (clinically & radiographically).
- Large diameter suction tips are not recommended as they can spread the gel to other areas.
- A patient consent form is available for download at www.sdi.com.au

**PROCEDURE:**

- Take a pre-operative photograph to record tooth shade.
- Take pre-operative X-ray to verify the presence of an acceptable root canal treatment and the absence of apical pathology.
- Isolate tooth/teeth undergoing whitening using a rubber dam. Ensure the dam fits tightly at the cervical margin.
- Open lingual access to pulp chamber and remove all restorative material and ensure correct endodontic access has been achieved. Remove gutta-percha to 1-1.5 mm below the cemento-enamel junction (CEJ).
- Place a 1 mm thick protective barrier eg a conventional glass ionomer on the exposed gutta-percha.
- Take a periapical radiograph to check that the barrier has been well placed.
- Clean pulp chamber with water.
- Firmly** attach a mixing nozzle to the Pola Professional 37.5% syringe away from patient. Dispense a small amount of gel onto a mixing pad until a uniform gel is extruded.
- Gently extrude a small amount of Pola Professional 37.5% gel into the open pulp chamber.
- Leave gel in for 8 min (Optional: curing light may be used at this stage, please refer to light manufacturer's recommendations) NOTE: If using a curing light, two curing lights can be used (one on each side of tooth) to activate Pola Professional 37.5%.
11. Suction gel out.
12. Repeat steps 9 to 11 for up to 4 applications.
13. Coat the inner labial surface of access cavity with bleach gel, place dry cotton pellet. Seal with composite material (we recommend Glacier (SDI Ltd) in a contrasting shade) without using a bonding agent (this helps prevent excess pressure build up within tooth).
14. Verify occlusal clearance.
15. Recall patient 2-3 days later for assessment. Procedure can be repeated until desired shade is achieved - first place rubber dam, remove temporary cover and gel then repeat steps 13 and 14. Note: The longer the tooth has been discoloured, the longer it will take for the bleaching treatment to remove the discolouration.
16. Remove gel from pulp chamber and wash with water to remove residue. **Note: Gel is quite sticky, ensure chamber is free from all gel residue.**
17. Place a temporary seal, (we recommend a glass ionomer cement such as Riva Protect (SDI Ltd) in a pink shade).
18. Wait 10-14 days before completing final restoration.
19. Record final shade.

**STORAGE AND HANDLING:**

- Store in cool place (2°-25°C / 35°-77°F).
- Do not use after expiration date.

**PRECAUTIONS:**

- For professional in-office use only.
- Patient must wear protective eyewear.
- Operators must wear gloves, mask and protective eyewear.
- Do not use on pregnant or lactating women.
- Do not use on children under 14 years of age.
- Do not use the Gingival Barrier on any persons having known resin allergies.
- Do not use Pola Professional 37.5% on any persons having known peroxides allergies.
- Anyone with a history of chemical allergies is advised to undergo allergy checks before treatment.
- Single patient use only - each syringe must not be used on more than one patient.
- Do not allow Pola Professional 37.5% gel to come in contact with skin, eyes and soft tissues.
- Pola Professional 37.5% will not lighten any restorative materials.

**FIRST AID:**

- SKIN/TISSUES (contact):** Immediately wash thoroughly with water. Apply a slurry of sodium bicarbonate and water to the affected area. Seek medical attention if symptoms persist.
- EYES (contact):** Open the eye wide and thoroughly wash for 15 minutes with running water. Seek medical attention if symptoms persist.
- INGESTION:** Rinse mouth with water. Gargle with salt water and drink lots of milk or water. Seek medical attention if symptoms persist.
- INHALATION (powder):** No symptoms expected.

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