· Conscious dentistry ·

# VIDU

## FLOWABLE COMPOSITE

es · COMPOSITE FLUIDO  
it · COMPOSITO FLUIDO  
fr · COMPOSITE FLUIDE  
pt · COMPÓSITO FLUÍDO  
de · FLIESSFÄHIGES KOMPOSIT

“THE FUTURE  
IS IN YOUR HANDS”

# VIDU

FLOWABLE COMPOSITE VIDU

en

### INSTRUCTIONS FOR USE

FLOWABLE COMPOSITE VIDU is a light cure, radiopaque flowable composite material.

#### INDICATION:

Class III, IV and V restorations.

#### RECOMMENDED PROCEDURES:

1. Clean out and dry preparation. Choose shade. 2. Prepare the cavity according to the general rules for adhesive dentistry. Bevel enamel margins. Protect pulp in deep cavities, according to general rules of preparative dentistry. 3. Etch enamel margins with Etch for 20 to 35 seconds. 4. Rinse off etching gel completely. 5. Dry with oil-free air. Avoid dentin desiccation. 6. Apply the adhesive (Bonding) according to its instructions. 7. Select the syringe containing the desired shade of FLOWABLE COMPOSITE VIDU. Remove syringe cap and insert an angled luerlock syringe tip. Lock tip into place. 8. Pressing the syringe plunger dispenses FLOWABLE COMPOSITE VIDU. To ensure a complete cure do not place FLOWABLE COMPOSITE VIDU in layers exceeding 2mm. 9. Light cure each layer of FLOWABLE COMPOSITE VIDU for 30 to 40 seconds. 10. Repeat steps 8 and 9, until the cavity is filled. 11. Remove excess, shape and finish with diamond burs and flexible disks of increasing grits.

#### ESSENTIAL NOTES:

Light curing units should have 450nm at 500mW/cm2 output and should be checked regularly. Place the light as close as possible to the restorative material. Do not disturb the oxygen inhibited layer because it is needed for bonding to the next layer. The use of a matrix is advised. Follow the manufacturer's instructions for use of any materials which are used with FLOWABLE COMPOSITE VIDU.

#### SAFETY PRECAUTIONS:

Placement is contraindicated, if a dry working area or the described application technique is not possible. Do not use on the exposed pulp or if an allergy to any component is known. Side effects are not known to date. In singular cases hypersensitivity to components of the material cannot be excluded. Avoid contact of paste with skin. If accidental contact occurs, wash immediately with soap and water. In case of contact with eyes rinse with copious amounts of water and consult a physician immediately. Avoid materials containing eugenol, moisture and oily air; they can inhibit the polymerization at the contact area. The use of cationic mouth washes, plaque indicators and chlorohexidine may lead to discolorations.

#### SIDE EFFECTS:

No side effects are known to date. Hypersensitivity against components of the material cannot be excluded in every case. Avoid skin contact with the paste. In the event of inadvertent contact, wash immediately with soap and water. In the event of eye contact, immediately rinse with plenty of water and seek medical treatment. If you experience any undesired effects, even in case of doubt, we ask for notification in every case, describing the concomitants and symptoms as detailed as possible. Product should not be used during pregnancy.

#### PROCESSING INSTRUCTIONS:

Observe the processing instructions in the instructions for use included with the product and the precautionary measures under all circumstances. Should the patient still show an undesired reaction or if there are undesired effects, we ask for notification in every case, describing the concomitants as detailed as possible. Recommendation: Grinding outside the mouth only under an intact fume hood.

#### COMPOSITON:

Dental glass in an optimized BIS-GMA (oligomere Dimethacrylates) resin matrix, containing 63% inorganic fillers. Product contains titanium dioxide less than 0,1 % (G/G)

#### STORAGE AND SHELF LIFE

Store dry at <25°C (78°F).  
Do not use after expiry date.

#### PRESENTATION:

- Syringe x 1ml - Shade: A1 [REF. V192101], A2 [REF. V192102], A3 [REF. V192103]  
A3.5 [REF. V192104], B1 [REF. V192105] and B2 [REF. V192106].

Keep away from children!  
For professional dental use only!

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