· Conscious dentistry ·

VIDU

GLASS-IONOMER CAVITY LINER

es · IONOMERO DE VIDRIO FONDO CAVITARIO
it · VETROIONOMERO SOTTO FONDO
fr · VERRRE IONOMERE FOND DE CAVITÉ
pt · IONÓMERO DE VIDRO FORRADOR
de · GLAS IONOMER UNTERFÜLLUNGSMATERIAL

“THE FUTURE
IS IN YOUR HANDS”

VIDU

GLASS-IONOMER CAVITY LINER

en

INSTRUCTIONS FOR USE

Product description

Light-cured, radiopaque, single- component underfilling material.

Indication

- Underfilling material for composite restorations
- For blocking out undercuts

Basic safety information

- Only for use by dental care professionals
- Store in places inaccessible to children

Application

1. Before preparation, clean the teeth and remove all cleaning residues. If the dentine over the pulp is thinner than 0.5 mm (pulp showing through), the pulp must be capped (with a calcium hydroxide preparation) to neutralise any potential clinically relevant diffusion of synthetic components.
2. Prepare cavities according to the general practice for adhesive techniques. Protect pulp in deep cavities using appropriate materials.
3. Etch cavities using your preferred technique (e.g. total etch). Follow the manufacturer's instructions for the materials used for this.
4. Apply the bonding agent.
5. Put the required amount of material on a mixing pad shortly before application. Protect the material from daylight and operating light to prevent the material curing prematurely. Apply a thin layer of the underfilling material using a small, round condenser. The layer should not exceed a thickness of 1 mm. If a thicker layer is required, apply the material layer by layer; allow each layer to cure before applying the next. The material may only be applied on the dentine and not on the enamel. Remove excess material and use a halogen light (at 450 nm with an intensity of at least 400 mW/cm²) to cure the material.
6. Bring the light source as close as possible to the underfilling material while doing so. Allow the underfilling to cure for at least 30 seconds with the light source approx. 2 mm from its surface. Never allow the material to cure for less than 20 seconds. Any stubborn residues on the enamel should be removed using a rotary instrument.
7. Apply the filling material (composite) according to the manufacturer's instructions.

Addend

The liner is so strong after curing that a probe cannot be used to remove excess material from cavity linings. This would loosen the entire lining. Instead, use a rotary instrument to carefully remove the excess. The material should be stored at room temperature, but not exceeding 22°C. Refrigerating the material (4-8°C) will increase its shelf-life. Allow the material to resume room temperature before use.

Side effects

Undesirable side-effects to this medical-product are not to expect, if working corresponding to the instructions of use. In case of side-effects - even when doubtful - please inform us with exact description circumstances and symptoms.

Contraindications

Do not use on patients with known allergy to the product's ingredients. Interaction with other materials. The liner should not be used in conjunction with materials containing eugenol, as eugenol impedes its setting reaction.

Instruction references

Take care of all working- and caution advices. In case of unexpected side effects or other reactions (even if work was made according to the Instructions of use) please inform us, if possible, with exact description of the accompanying circumstances and symptoms.

NOTE

Applicator-tools (tips) are only determinate for single-use at only one patient.

Composition

Methacrylate and glass-ionomer powder

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