COMPOSICIÓN:

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
COMPOSITION:

100 g powder containing 93 g zinc oxide.
100 g liquid containing 61 g phosphoric acid.

FORTEX cement is a product used in dentistry for fixing crowns, bridges, fillings, etc. and also for prosthetic works in orthodontic work-shops.

This cement consists of two components:

- A powder based on zinc oxide.

- A liquid based on phosphoric acids.

The reaction of these two components when mixed is to neutralize themselves and to form a rock hard mass of zinc oxyphosphate.

Instruments required for the mixing:

In order to mix the components of **FORTEX cement** a flat, clean and dry surface is needed, for example a glass tile, and a suitable stainless steel spatula. Mixing proportions: 0,5 g of **FORTEX** powder to 5-6 drops of **FORTEX** liquid. (1 part weight of the liquid to 1,5 part weight of the powder).

Mixing process:

The mixing time should take from 45" to 1 minute.

Take the appropriate quantity of **FORTEX powder** and liquid for the pieces to be bonded. Mix a small quantity of powder with the liquid using quick, continuous and circular movements (only use small quantities of powder) until obtaining a creamy and even mass which will drip.

The paste is now ready to be used for fillings or for coating the parts to be bonded and keeping them in place until the **FORTEX cement** has completely set and hardened.

The filling consistency is reached by increasing slightly the proportion of powder and when, on lifting the spatula, the extended point folds over to form a small hook and does not sink back into the mixture (Fig. 1). In this case, the setting period is 3-4 minutes.

The Cementing consistency is reached with the normal mix proportions and when, on lifting the spatula, the point so formed slowly sinks back into mixture. In this case the setting period is 4-6 minutes (Figs. 2, 3, 4).

Setting of the FORTEX cement:

The setting time of the **FORTEX cement**, after which it cannot be worked, is of some two and a half to three minutes, provided that specific working conditions are respected:

- A mixture proportional to 0,5 grams of powder and 6 drops of liquid.

- An ambient temperature around 15°C - 20°C.

- The temperature of the surface to be worked and the instruments must not greatly differ from the above temperatures.

Recommendations:

If you want to reduce or increase the setting time, this can be done by making the paste thinner or thicker: By adding a larger quantity of powder than that recommended the setting process is speeded up. By reducing it, the setting process is delayed. Do not use too much liquid, it could be that the product will not set or it could irritate sensitive tissues such as dentalpulp. Remember that the surface to be worked must be completely clean and dry in order to facilitate the total bonding of the pieces used.

When working on sensitive surfaces, these must be adequately protected.

For lining silicate and composite fillings the silicate or composite is pressed on the **FORTEX** cement mixture when still soft.

Preservation:

The **FORTEX cement** must be protected from the air and humidity and it is therefore recommended to close the bottle immediately after use.

FORTEX TRIPLE: 3 x 100 g of powder.
3 x 30 ml of liquid.

FORTEX CLÍNICA: 100 g of powder.
30 ml of liquid.

FORTEX POLVO: 100 g of powder.

FORTEX LÍQUIDO: 30 ml of liquid.

FACIDEN, S.L. accepts no responsibility for preparation and use in ways other those described in this leaflet, nor does it accept responsibility for packages sold without the pertinent seal.

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