Amalgama
non gamma 2

ventura high alloy

Amalgam
non gamma 2

CE
2460

Rev 09

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1. DESCRIPCIÓN

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EN INSTRUCTIONS FOR USE

1. DESCRIPTION

Ventura high alloy is a dispersed phase amalgam alloy: a mixture of lathe-cut particles and microfine eutectic spherical particles. The handling characteristics are comparable to a conventional amalgam: sufficient condensation pressure should be applied.

Ventura high alloy is free from the corrosion-prone and weak gamma-2 phase. Durable restorations with a lasting lustre can be made.

Alloy composition: Silver (Ag) : 69.2% (m/m),
Tin (Sn) : 18.6% (m/m),
Copper (Cu) : 11.9% (m/m),
Zinc : 0.3% (m/m).

Ventura high alloy is in full compliance with the latest version of the ISO24234. The mercury contained in the pre-dosed capsules is in full compliance with ISO 24234:2015.

2. INDICATION

Stress-bearing fillings in posterior teeth (Class I and II).

3. CONTRAINDICATIONS

- Use of amalgam is not recommended when a physician has deemed it unsuitable for the patient.
- Use of amalgam is not recommended for the treatment of deciduous teeth, children under 15 years and pregnant or breastfeeding women, except when strictly deemed necessary by the practitioner on the ground of specific medical needs of the patient (from 1 July 2018).
- Hypersensitivity (allergic reaction) to amalgam components. In such cases the use of suitable alternative materials must be considered.

4. WARNINGS

- After placement or removal of dental amalgam restorations increased mercury concentration in blood and urine has been observed. According to available scientific knowledge this increase has not been associated with any adverse health effects.
- If placed in close contact with other metal restorations, galvanic effects may occur. In most cases they will be of short lasting duration. If the effect persists, the user should consider replacement of the dental amalgam filling with another material.
- In individual cases, local mucosal reactions (lichenoid) have been observed. Such local reactions may be of an irritative (mechanical, chemical, electrochemical) or allergic nature. In case of allergy to components of dental amalgam the use of suitable alternative materials must be considered.
- Unnecessary exposure to mercury vapour or dental amalgam particles during handling, placement or removal of dental amalgam should be avoided. Placement and removal of dental amalgam fillings should be performed with appropriate water spray and vacuum suction.
- Instructions and regulations for storage and disposal of dental amalgam waste must be observed.

5. DIRECTIONS FOR USE

Mixing
Ventura high alloy is available as pre-dosed capsules. It should be noted, that the working- and hardening time would be somewhat different, dependent on the type of mixer used: a high-energy mixer gives a faster hardening. The mixed amalgam should be a plastic, coherent mass with a shiny surface. Mixing ratio for **ventura high alloy** is: 10 parts alloy / 11 parts mercury by weight.

Ventura high alloy in all cases 5 – 7 seconds in High-energy mixer (>4200 rpm) such as **ventura mix 2**.

	Capsule	Alloy	Hg
VENTURA HIGH ALLOY CAP 1	PINK SKY BLUE	400 mg	435 mg
VENTURA HIGH ALLOY CAP 2	PINK BLUE	600 mg	653mg
VENTURA HIGH ALLOY CAP 3	PINK WHITE	800 mg	865 mg

Nota:

1. For safety reasons it is best to reclose the empty capsule after taking out the mixed amalgam, before cleaning or disposal.
2. The pre-dosed capsules (Type 1) are self-activating: contact between alloy and mercury occurs automatically when the capsule is subjected to the action of the mixing machine. Length: 32,2 mm, max. exterior diameter: 13,7 mm.
3. Storage at temperatures no higher than 25°C is recommended.
4. Dispose of cotton wool and used capsules in accordance with national regulations.

6. INSERTION AND CONDENSATION

After mixing, the amalgam should immediately be inserted into the cavity using an amalgam carrier. The filling should be built up in layers under adequate condensation pressure using a plugger appropriate for the size of the cavity. Do not use an ultrasonic condenser. Carefully remove excess mercury.

7. FINISHING AND POLISHING

Approx. 8-10 minutes after mixing, the filling can be finished. Polishing should be done in a subsequent session.

THIS ALLOY CONTAINS ZINC; THE AMALGAM MADE THEREFROM MAY SHOW EXCESSIVE EXPANSION IF MOISTURE IS INTRODUCED DURING MIXING OR COMPACTING.

8. LABELLING

This information is in compliance with ISO 24234:2015 and also with European legislation.

Danger

9. RECOMMENDATIONS FOR SAFETY PRECAUTIONS

Dental clinical personnel, when using silver amalgam, should exercise proper procedures to avoid personal contact with mercury. Precautions should be taken to avoid exposure to mercury vapour in the dental environment. Waste amalgam should be collected for disposal in a manner, which will protect both those who handle the waste and the environment.

1. Treat mercury and amalgam with care. Be MERCURY-AWARE.

2. Ideally a high-energy mixer should be equipped with a protective cap to cover the capsule whilst mixing.
3. Small droplets of spilled mercury can be collected with the aid of a 'mercury-collecting forspes' or a small amount of freshly mixed amalgam that will easily absorb liquid mercury. NEVER use a vacuum cleaner!
4. Placement and removal of amalgam fillings must be done under water-cooling and with effective vacuum suction. The use of eye protection and mouth mask is advisable.
5. The wastewater drain of the spittoon should be equipped with an amalgam separator. Local regulations on type and installation of such a separator and handling/recycling of the amalgam waste should be taken into consideration.
6. Excess mixed amalgam should be kept in a container, designed for the purpose, or at least in a closed container filled with photographic fixer.

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