• Conscious dentistry •

TULIP 100A

en • Dental curing light LED

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English (Pág 4-19)

1. **USE AND LIABILITY**
2. **PART NAME**
3. **INSTALLATION AND USE**
4. **GENERAL USAGE**
5. **VIGILANCE**
6. **SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**
7. **MANUFACTURER'S RESPONSIBILITY**

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1 • USE AND LIABILITY - RISK ASSOCIATED WITH THE PRODUCT

Prior to installation and start-up of the unit, please read these instructions carefully!

- Please strictly follow the instructions, we are not responsible for the consequences caused by improper use.
- The light curing two luminous time interval is not less than 60 seconds.
- The light guide rod is fragile glass products, please use carefully, when the light guide rod has obvious damage, cracks will affect the curing effect, please replace the new light guide.
- Disinfection method of light guide rod: Soak with 75% alcohol for a few minutes or heat with 121°C high temperature wet heat disinfection
- After full charge, please remove the charger power plug from the power socket.
- When used for a period of time (more than 80 times), the battery voltage is reduced, which will affect the curing effect, and should be charged in time. Due to the use of high-performance the lithium-ion battery can be charged at any time without affecting the service life of the battery.
- Battery electrodes at both ends do not touch the metal, so as to avoid short circuit damage to the battery.
- Please regularly wipe the contact point of the battery and the charging needle on the charging seat with a cotton swab dipped in a little anhydrous alcohol, so as to avoid the battery caused by poor contact or the occurrence of charging dissatisfaction.
- To prevent blue, ultraviolet and thermal radiation, do not look directly at the light source, and wear safety goggles when operating.
- Please do not use the power supply other than the original charger, otherwise it may cause serious damage to the equipment and personal injury.

- Unauthorized disassembly and modification will lead to the expiration of the warranty, the company is no longer responsible for the product.
- In addition to replacing the light guide rod and rechargeable battery, it needs to be maintained and repaired by professionals.
- Do not use this equipment in strong magnetic field environment.
- Do not use this equipment in lightning weather.

BATTERY CHARGING

When the battery is found to be insufficient, it must be charged or replaced with a fully charged battery.

Charging is required in any of the following cases:

- The host has an alarm
- When buying a new battery
- Not using the battery for a long time

NOTES:

- The host is equipped with undervoltage detection function, when the battery voltage is lower than the set value, the host will send an alarm.
- Please charge in time when automatic shutdown occurs to avoid battery failure due to excessive discharge.
- When purchasing a new battery, due to insufficient charge from the factory or storage for a long time, the battery voltage has a low voltage value, so charge it in time.

- Generally charged once, not used for a long time or has been used several times, should be recharged in about a month.

PRECAUTIONS FOR USING THE BATTERY HOST:

- During use, keep away from heat source and high pressure to avoid children playing with the battery.
- This battery can only be charged using the matching charger, do not charge the battery in the charger for more than 24 hours.
- Do not short circuit the positive and negative electrodes of the battery, and do not let the battery be damp to avoid danger.
- When the battery is not used for a long time, it should be removed from the device and stored in a low temperature, cool and dry environment.
- Dispose of waste batteries safely and properly, and do not throw them into rubbish bins, fires or water.

SAFETY CLASSIFICATION OF EQUIPMENT

- Classified by anti-shock type: Class II equipment with internal power supply
- Classification according to the degree of shock protection: BF type application part
- Classified according to the degree of protection against liquid intake: ordinary equipment (IPX0), not waterproof
- Operation mode classification: intermittent operation
- Degree of safety in the case of burning anesthetic gas: can not be easily mixed with air
- Equipment used in the presence of anesthetic gas or flammable anesthetic gas mixed with radon or nitrous oxide
- Function and intended use: Suitable for dentistry, it has the function of accelerating the curing of dental repair materials.

CAUTIONS FOR USERS AND PATIENTS

- Direct or indirect exposure of the eyes should be avoided. Prolonged exposure of the lamps is unpleasant for the eyes and can cause damage. Therefore, the use of the included anti-glare cones is recommended. Individuals who are generally sensitive to light, those taking photosensitizing medications or those who have undergone eye surgery, or persons working with or in close proximity to the device for long periods of time should not be exposed to light from this device and should wear protective goggles (orange) that absorb light below 515 nm.
- CAUTION. Restrict exposure to the area of the oral cavity in which clinical treatment is intended. Protect patient and user from reflection and intensive scattered light by taking the appropriate precautions, e.g., glare shields, goggles, or coverings
- As with all high-performance lamps, the high light intensity causes some thermal development. Prolonged exposure of areas close to pulp and soft tissue can cause irreversible or reversible damage. Therefore, the recommended polymerization times should be observed and uninterrupted polymerization times of more than 20 seconds on the same tooth surface should be avoided, as well as direct contact with the gingiva, dental mucous membrane or skin. If necessary, polymerize indirect restorations in intermittent intervals of 20 seconds at a time or use external air-gun cooling.
- If it is suspected that the unit is no longer safe to operate, it must be disconnected from the mains to prevent accidental operation. This could occur, for example, if the device is visibly damaged or no longer functions properly. Complete interruption of power is only ensured by unplugging the device from the mains.

CONTRAINDICATIONS

- Materials whose polymerization is activated outside the wavelength range of 385-515 nm (no material known to date). If you are

unsure about certain products, consult the manufacturer of the respective material.



- Do not charge or use the device near flammable or combustible substances.
- Portable or mobile high-frequency communication devices may interfere with medical equipment. The use of cell phones during operation of the lamp is not permitted.

2 · PART NAME

1. Handle
2. Head
3. Light shelter
4. Charger
5. Power supply

1. Open/Stop
2. Model screen
3. Model/time setting
4. Battery

1. Charge seat

2. ILgh intensive tester
3. Charge screen

3 · INSTALLATION AND USE

INTENDED USE

Tulip 100A adopts the principle of ray radiation to solidify the light-sensitive resin by firing it in a short time. This product is used in dentistry. It has the function of accelerating the polymerization of dental restorative materials (adhesives, composites, ionomers, cements) and solidifying tooth whitening material. The beam emerging from the device covers the light wavelength range of 430 to 480 nm. The intended place of application is in the dental office, doctor's office or hospital.

WARNING: Only use the handle and head supplied with this product.

1. As shown in Figure 1 below, the machine head can be pulled out and inserted into the handle along the direction of the arrow; Insert the head of the machine requires a little force, when the head is completely stuck into the handle, there will be a crisp "click" sound; Pull out the head with a little force when you need to remove it.
2. As shown in Figure 2 below, it is the state diagram of the connection between the handle and the machine head. The machine head can rotate 360 degrees, and the appropriate treatment Angle can be selected according to the use habits and the actual situation.


Fig 1.


Fig 2.

3. Installation and removal of light shield

NOTE: Before installing the light shield, a disposable isolation sleeve should be used to cover the head of the handle to avoid contact between the head or other components and the skin or oral mucosa of the patient. Disposable condoms are not reusable and should be disposed of in an appropriate manner according to local laws and regulations.

As shown in Figure 3 below, the gap of the light shield can be installed and removed along the direction of the arrow, and the gap of the light shield can be stuck in the appropriate position of the machine head according to the first step of the figure. Then according to the second step to the direction of the handle sliding light block board, until the light block board is tight and not loose.

As shown in Figure 4 below, it is the finish of the installed light shield.


Fig 3.


Fig 4.

4. The charging base is connected to the power supply

WARNING: Use only the charging base and power supply provided with this product.

- Connect the power supply plug to the jack on the back of the charging base as shown in Figure 5 below.
- Plug the other end of the power adapter into AC100,240V socket

5. Handle is connected to the charging base

WARNING: Only use the handle and charging base provided with this product.


Fig 5.

- As shown in Figure 6 below, when the battery of the handle is dead or needs to be charged, put the battery compartment of the handle into the one connected to the power supply
- The handle charging jack of the charging base can charge the handle.
- When charging, there will be an icon in the display window of the charging holder (refer to the system icon for details)

NOTE: There is a battery low power detection function inside the handle. When the battery low power is detected, the power symbol on the screen of the handle will display "0" power and flash. Please charge in time.

6. Replacement of batteries

WARNING: Only use the batteries provided by the after-sales service or the dealer, do not use other brands of batteries, so as not to cause damage to the light curing machine or other hazards; Please confirm whether the battery needs to be replaced according to the actual situation. It is strictly forbidden to remove the battery at will.


Fig 7.


Fig 6.

- If you need to replace the cone battery, first rotate the bottom cover of the handle end in a counterclockwise direction, then remove the cone battery and replace it with a new battery, and finally screw the bottom cover into the handle in a clockwise direction.

7. System display

Display	Model	Time setting (seconds)	Light intensive (mw/cm2)
 10s M1	M1 standard	5,10,15,20	1000-1300
 05s M2	M2 strong	3,4,5	1800-2100
 03s M3	M3 super strong	1,2,3	2700-3100
 15s M4	M4 orthodontics	3 secs*5; 3 secs*10	2700-3100
 30s M5	M5 decay detective	30,60	/
 10s M6	M6 pulse	5,10,15,20	1000-1300
 20s M7	M7 increasing	5,10,15,20	1000-1300

8. Charge base display

Charge base display	Name	Usage
	Charging	When charging the handle, the charge base will display it
	Light intensive	When measuring light intensive, the charge base will display it

9. How to select models and time.

Key	Name	Usage
	Power on/start/stop	Press one time to power on; then press one time to start; then start one time to stop
	Model/time	Use this button to setting model and time. Short press to select models, long press to select time.

10. Measuring light intensity

The charging base has a built-in light measuring core, which can quickly measure the light intensity when the charging base is connected to the power supply. Select the corresponding model of the handle, and place the light source of the head in the center of the metering area.

NOTE: Placing the head light source at an oblique Angle to measure the light intensity will lead to a significant deviation in the measured light intensity value.

The effective measured light source diameter of the built-in light metering core is 10mm, if the measured light source diameter is not 10

mm(as is common 8mm), resulting in significant differences in the measured light intensity values. At the same time, there are differences in the standard of measuring the light intensity value of different manufacturers, it is recommended to only use the matching light curing machine. If the measurement is not the matching light curing machine, the measured light intensity value may be significantly different.

4 • GENERAL USAGE

When using, first put the disposable isolation sleeve into the head of the machine, and the light shelter into the grip, and then the light source is aligned to the position that needs to be cured (the distance between the smooth surface and the cured surface is recommended to be small 1000cm), press the button , the handle will emit a sound of "beep", the head of the machine will emit blue light and work according to the selected working mode, the handle will display the countdown time from the selected timing time, and the work will end when the timing reaches 0, and the handle will display the return to the set timing time.

In the normal use of the time before the end of the timing, press the button  can end the work at any time to turn off the light source.

1. Troubleshooting

If you meet any problem during usage, please contact our local distributor or us.

2. Product technical parameters

HANDLE AND NOSE

Overall assembly size: 22.5mm X 22.5mm X 205mm

Overall assembly weight: 110 g

BATTERY

Battery type by cell: Rechargeable batteries

Battery model: 18500 3.7V

Battery capacity: 2000mAh

Battery protection: battery with over voltage, over current and short circuit protection

Power adapter (charger)

Input: 100-240V 50/60Hz 0.4A Max.

Fuse breaker: T1A250V

Output: 5V $\pm$ 5% 1 A

THE PERFORMANCE OF LED LIGHTS

Wavelength: 385nm~515nm

The typical wavelength peaks of the LED lamp of the photosolidification machine are: 465 nm and 400 nm

Category: Class I

Emission limit (AEL) : 3.0×10^{-3} W/cm²

Inspection method: when the LED light is used correctly, the light is in good condition.

Surface diameter of light emitting element: 10 mm

Optical effective area: 78mm²

standard requirements 385nm~ 515nm (blue light), not less than 250mW/cm²

EQUIPMENT SAFETY CLASSIFICATION

Classification by type of electric shock: Class II equipment with internal power supply

Classification by degree of electric shock: Type B application part

Classification according to the degree of protection for the liquid: Ordinary equipment (I PXO), not waterproof

Classification according to the operation mode: Intermittent operation

ENVIRONMENT OF USE

Ambient temperature: - 20°C~+20°C

Relative humidity:10%~93%

Air pressure: 70KPa~106KPa

5 • VIGILANCE

If any serious incident has occurred in connection with the device, inform the manufacturer and the competent authority of the Member

State in which the user and/or the patient is established.

6 • SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

TULIP 100A is safe and performs as intended when used in accordance with the manufacturer's instructions for use. The user shall be responsible for determining if the product is suitable for his/her application.

7 • MANUFACTURER'S RESPONSABILITY

Our products have been developed for professional dental use. Since the application of our products is beyond our control, the user is fully responsible for it. Of course, we guarantee the quality of our products according to the applied standards.

	Caution		European Union certification
	Manufacturer		European Union Representative
	Sterilization Temperature		Heat Sterilized
	Serial Number		Manufacturing Date

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Being Foshan Medical Equipment Co., Ltd
Shanghai International Holding Corp. GmbH