VOCO**BrackFix®**

Light-curing bracket bonding system

Instructions for use**Product description:**

The **BrackFix** bracket bonding system consists of the light-curing fluorescent **BrackFix Adhesive** and a primer which is available as an etch & rinse or self-etch primer. **BrackFix** is used for the direct bonding of orthodontic metal and ceramic brackets. The individual components can only be combined with light-curing materials.

Indications:

Direct bonding of orthodontic metal and ceramic brackets

Application:**A. Preparation of the teeth**

- When using **BrackFix Primer** (etch & rinse) or **BrackFix Primer SE** (self-etch), examine the enamel surface of the tooth areas to be adhered for structural damage (cracks, lesions, and the like). Damage of this type can cause enamel to be pulled off when the brackets are removed. When using **BrackFix Primer SE**, attention should also be paid to fluorosis and discolourations. Fluorosis can lead to impairment of the adhesive strength. The patient should be accordingly informed prior to the start of treatment. Performing the bonding on ceramics is not recommended as otherwise damage to the ceramics cannot be excluded.
- Isolate teeth.
- Clean teeth with pumice powder or polishing paste. Rinse with water.
- Carefully dry with dry, oil-free air.

B. Conditioning the teeth (when BrackFix Primer is used)

- Before using **BrackFix Primer** (option 1), etch the enamel surfaces with a phosphoric acid gel. In doing so, proceed according to the instructions for use of the etching gel.
 - Rinse with water.
 - Carefully dry with dry, oil-free air.
 - Avoid recontamination with saliva and water.
- When using **BrackFix Primer SE**, follow option 2 (see below).

C. Priming the teeth**Option 1: Applying the etch & rinse primer (BrackFix Primer)**

Measure out **BrackFix Primer** into a mixing palette and then use a suitable applicator to apply as a thin homogeneous film on the conditioned tooth surface. Note: The primer is light-curing and therefore excessively intense ambient light exposure should be avoided. Dim operating lights during the application.

Option 2: Applying the self-etch primer (BrackFix Primer SE)

- Detach a blister of *SingleDose* at the perforation and turn the printed side up. Hold the blister of *SingleDose* between thumb and forefinger and, by pressing on the area marked "press here", allow the liquid contained in the blister to flow into the mixing and dispensing chamber. Apply the enclosed *Single Tim* applicator to the centre of the coloured circle in order to pierce through the film of the mixing and dispensing chamber. Expand the hole to its maximum size using a circular motion. By stirring thoroughly with the applicator, create a homogeneous, streak-free mixture of the two liquids.
- Apply **BrackFix Primer SE** in a thin, homogeneous layer on the tooth surface and rub in for 3 - 5 s.
- Remove more material from the *SingleDose* and perform the same procedure on the next tooth.
- When all teeth in a dental arch have been coated with the Primer SE, dry the Primer SE on each tooth with a soft, oil- and water-free air stream in a mesial-to-distal direction. Proceed quickly with the application of **BrackFix Adhesive** (see below). If the process is delayed, apply another layer of Primer SE and proceed accordingly.

D. Application of adhesive

The **BrackFix Adhesive** is applied after the primer is applied to the (preconditioned) enamel surfaces. To do this, apply a sufficient amount of **BrackFix Adhesive** to the base surface of the bracket. Remove any remaining material protruding from the syringe opening and close the syringe with the lid.

E. Application and polymerisation

- After applying **BrackFix Adhesive**, quickly position and bond the bracket on the surface of the tooth.
- Remove excess material at the edge of the bracket carefully using an appropriate instrument (for example, a probe).
- For the curing, use a commercial light-curing unit with a brightness of at least 1000 mW/cm². In the case of metal brackets, position the light guide as near as possible to the interproximal bracket surface and expose 10 s mesially and 10 s distally per bracket. In the case of ceramic brackets, position the light guide perpendicular to the bracket surface and expose as near as possible through the bracket for 20 s.
- Dental archwire can then be used.

F. Removal

Follow manufacturer's instructions for removing the brackets.

G. Final polishing

After removing the brackets, any remaining adhesive residues on the surface of the tooth can be made visible using a UV lamp (adhesive detectable at: 360 - 370 nm, min. 20 mW/cm²) and removed using a carbide tipped finishing burr (approx. 4000 rpm), for example. It is recommended that polishing and fluoridation be subsequently performed.

Information, precautionary measures:

- Phenolic substances (e.g. temporary materials containing eugenol) impair polymerisation and therefore must not come into contact with **BrackFix**.
- In case of eye contact with **BrackFix**, rinse eyes immediately with plenty of clean water and seek medical advice.
- The **BrackFix** adhesive system is not recommended for use with plastic brackets (polycarbonate).
- BrackFix Adhesive** and **BrackFix Primer** contain methacrylates, amines, and butylhydroxytoluol. Do not use on patients with a known hypersensitivity (allergy) to these ingredients.

- BrackFix Primer SE** contains organic acids, methacrylate, amines, and butylhydroxytoluol. Do not use on patients with a known hypersensitivity (allergy) to these ingredients.
- UV light can cause damage. Patients and dental personnel should always take note of the corresponding safety measures. This includes wearing appropriate protective UV goggles, among other things. Do not look into the light source and do not aim the light directly into the eyes. Only use the light at the treatment site.

BrackFix Primer SE

- The blister must not be bent or folded prior to activation. Should the blister exhibit apparent signs of damage such as folds or kinks, it must no longer be used.
- The film of the blister and the label must not be damaged by the fingernails either on its top or underside.
- No tools must be used in the activation. Only the thumb and forefinger may be used for activation, whereby the thumb must be placed exactly on the marking of the label and the forefinger directly on the opposing underside.
- The activation and the subsequent piercing with the *Single Tim* should not be performed near the eyes.
- Never press on both chambers at the same time or push the liquid to and fro between the depressions.
- Contact with **BrackFix Primer SE** may cause bleaching of the oral mucosa as a result of the coagulation of proteins. This is a reversible effect which subsides within a few days.
- Once it is opened, **BrackFix Primer SE** is light-curing and further processing must therefore take place immediately. Avoid intense lighting during application.

Storage instructions and application method:

- Close the **BrackFix Primer** bottle tightly after use. Store upright to ensure backflow of the liquid.
- Store **BrackFix Adhesive** and **BrackFix Primer** at temperatures of 4 °C to 28 °C. Bring to room temperature prior to use.
- Store **BrackFix Primer SE** at temperatures of 4 °C to 23 °C. Bring to room temperature prior to use.
- Do not expose the material to direct sunlight or the operating light.
- Do not use the product after its expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1207 E1 0917 99 © by VOCO
VOCO GmbH