

# IYSERPEN PLUS ASER TECHNOLOGY



#### TECHNICAL DATA

## CODE: PADLHYSD001.1 **POWER SUPPLY:** 100-240 V ~ POWER FREQUENCY 50-60 Hz 1.5A @ 115Vac, 0.75A @ 230Vac MEDICAL CLASS / LASER CLASS: **INSULATION CLASS** II - Applied Part Type BF LASER SOURCE: 3W DIMENSIONS / WEIGHT: 167 x 18 mm (L x Ø) / 122 grams WAVELENGTH:

### **Accessories included**

- Base
- Handful
- USB power cable
- 10 Interchangeable Tips
- Large areas terminal
- Protective goggles
- Clinical protocols
- **User Manual**

Wireless semiconductor laser designed and certified for the dental clinic, preset with seven protocols which generally cover more than 80% of the daily activity:

810nm


**OPERATIVE MODE:** Selectable: Continuous -Timed

Biostimulation, Dental Pockets Decontamination, Desensitization, Whitening, Endodontics, Peri-implantitis, Herpes-Aphthae.

Hyserpen Plus represents the excellent compromise between usability , thanks to the display icon graphic, on the charging base, and the high level of technology, ensured by

multiple uses allowed by the device. Thanks to its compact size and with a weight of only 122 gr, Hyserpen Plus allows you to operate adequately on each patient and in many working conditions. Moreover you have the possibility to manage its even through a simple external power bank. This is the right solution for who wants to increase the productivity, rising up the number of devices available in the clinic, or for the professional that is approaching for the first time the use of laser technology

### Preset Protocols





SWISS & WEGMAN Srl