

NEUBER

Dental Clinic Cart with Glass Countertop

Clinical Furniture · Dental Equipment

Ref: UM-540



PRODUCT DESCRIPTION

The NEUBER UM-540 is a professional dental clinic cart designed for demanding clinical environments. It features a tempered glass countertop and a stainless steel structure, combining functionality with ease of cleaning and disinfection. The cart includes **5 units** (1 closing panel + 4 storage drawers): the top unit acts as a **closing/locking panel** with a security lock, concealing and securing the 4 functional storage drawers below. All surfaces are resistant to scratches, water and corrosion.

KEY FEATURES

■ 5 Units (1 closing panel + 4 storage drawers)	The top unit is a locking closure panel — NOT a storage drawer. It covers and secures the 4 functional drawers underneath.
■ Tempered Glass Countertop	Hygienic, scratch-resistant work surface, easy to disinfect with clinical-grade products.
■ Stainless Steel Structure	Corrosion-resistant body and guide rails ensure long service life in sterilization-intensive environments.
■ Security Lock	The closing panel is fitted with a key lock to restrict unauthorized access to the storage drawers.
■ High Cleanability	All materials — glass, stainless steel, aluminium handle and plastic tray — selected for easy cleaning and disinfection.
■ Compact & Mobile	500 x 500 mm footprint with castors for easy positioning alongside the dental chair or treatment unit.

TECHNICAL SPECIFICATIONS

PARAMETER	VALUE	NOTES
Dimensions (W x D x H)	500 x 500 x 840 mm	Compact dental cart format
Net Weight	32 kg	
Drawer Units	5	1 closing panel + 4 storage drawers
Countertop	Tempered glass	Scratch & chemical resistant
Body & Guide Rails	Stainless steel	Corrosion resistant
Handle	Aluminium alloy	
Tray	Plastic	Removable for cleaning
Color	White	
Security Lock	Yes — top closing panel	Key included
Water / Scratch Resistance	Yes	Splash & fluid resistant surfaces
Corrosion Resistance	Yes	Suitable for sterilized environments
Manufacturer Ref.	UM-540	Brand: NEUBER

■ **Important — Drawer Configuration:** The UM-540 has 5 drawer-height units. The **topmost unit is a closing/locking panel** — it does not open as a storage drawer. Its function is to cover and lock the **4 usable storage drawers** beneath it.