

## CLEANMED WIPES



- 200 wipes Dispenser



- 200 wipes Refill bag

Product Name	Code
Cleanmed Wipes Dispenser	7.004-D
Cleanmed Wipes Refill	7.004-R



### EFFECTIVE SPECTRUM

BACTERICIDAL	Contact Time
MRSA	1 min
E-Coli	1 min
Pseudomonas aeruginosa	1 min
Enterococcus Hirae	1 min
Staphylococcus Aureus	1 min
Salmonella Typhimurium	1 min
VRE	1 min
MYCOBACTERICIDAL	
Mycobacterium Avium	1 min
TUBERCOLICIDAL	
Mycobacterium Terrae	1 min
YEASTICIDAL	
Candida Albicans	1 min
VIRUCIDAL	
HIV1	1 min
H1N1 Influenza A	1 min
Hepatitis C Virus	1 min
Humana Influenza B Virus	1 min
Adenovirus	1 min
Norovirus	1 min
Leukemia Virus (HTLV)	1 min
Hepatitis Delta Virus (HDV)	1 min
Coronavirus (SARS, MERS, ...)	1 min
Filoviridae (including Ebola)	1 min
Rubella Virus	1 min
Measles Virus	1 min
Herpesviridae	1 min
Flavivirus	1 min
Rabies Virus	1 min
Poxviridae	1 min
Hepatitis B Virus	1 min
FUNGICIDAL	
Aspergillus brasiliensis	5 min

### PRODUCT DESCRIPTION

Alcohol based poly-sealed wipes, soaked with special formula for surfaces of non-invasive medical devices.

#### Main features:

- Alcohol based disinfectant (63% w/w / 70.5% v/v Ethanol);
- Benefit of the physical action ensuring the disinfectant makes full contact with the surface;
- **Wipes size: 200mm x 200mm**
- Virucidal against all enveloped Viruses including Coronavirus
- Effective against non-enveloped viruses Adenovirus and Norovirus
- Refill pack available
- Lemon fragrance

#### COMPOSITION

Aqua, Ethanol, Parfum, Linalool, Limonene.



#### SOLUTION

Ready to use



#### WARNINGS

- Do not use on alcohol sensitive surfaces



#### STANDARDS APPLIED:

- ✓ EN 13624: 2013
- ✓ EN 13727: 2012 + A1:2013
- ✓ EN 14348: 2005
- ✓ EN 14476: 2013 + A1:2015
- ✓ EN 14885: 2015