

CLEANMED WIPES AF



- 200 wipes AF Dispenser



- 200 wipes AF Refill bag

PRODUCT DESCRIPTION

Alcohol free poly-sealed wipes, soaked with special formula for surfaces of non-invasive medical devices.

Main features:

- Benefit of the physical action ensuring the disinfectant makes full contact with the surface
- **Wipes size: 200mm x 200mm**
- Virucidal against all enveloped Viruses including Coronavirus
- Compatible with materials and surfaces sensitive to alcohol
- Refill pack available
- Mint fragrance

COMPOSITION

Aqua, Benzyl-C12-18-Alkyldimethyl Chlorides, C12-14-Alkyl[(Ethylphenyl)Methyl] Dimethyl Chlorides, Didecyldimethylammonium Chloride, 2-Propanol (not an active, <1%), C9-11 Pareth-8, PPG-2 Methyl Ether, Polyhexamethylene biguanide hydrochloride, Cocamidopropyl Betaine, Tetrasodium EDTA, antifoam, preservative, aroma.



STANDARDS APPLIED:

- ✓ EN 13727: 2012 + A1:2013
- ✓ EN 14476: 2005
- ✓ EN 13624: 2003
- ✓ EN 14348: 2005
- ✓ EN 14885: 2015
- ✓ EN 13704: 2002



SOLUTION

Ready to use

Product Name	Code
Cleanmed Wipes AF Dispenser	7.004AF-D
Cleanmed Wipes AF Refill	7.004AF-R



EFFECTIVE SPECTRUM

BACTERICIDAL	Contact Time
MRSA	1 min
E-Coli	1 min
Pseudomonas aeruginosa	30 sec
Enterococcus Hirae	30 sec
Staphylococcus Aureus	30 sec
Salmonella Typhimurium	1 min
Legionella Pneumophila	5 min
VRE	1 min
MYCOBACTERICIDAL	
Mycobacterium Avium	2 min
SPORICIDAL	
Clostridium Difficile C-Diff	5 min
YEASTICIDAL	
Candida Albicans	30 sec
VIRUCIDAL	
HIV1	2 min
H1N1 Influenza A	2 min
Hepatitis C Virus	2 min
Humana Influenza B Virus	30 sec
Leukemia Virus (HTLV)	2 min
Hepatitis Delta Virus (HDV)	2 min
Coronavirus (SARS, MERS, ...)	2 min
Filoviridae (including Ebola)	2 min
Rubella Virus	2 min
Measles Virus	2 min
Herpesviridae	2 min
Flavivirus	2 min
Rabies Virus	2 min
Poxviridae	2 min
Hepatitis B Virus	2 min
FUNGICIDAL	
Aspergillus niger	30 sec