



**LEVEL 1
ISOLATION
GOWN
AAMI PB70**

TECHNICAL DATA SHEET

8060 | Isolation Gown

8060 Isolation Gown provides AAMI level I protection and is made from a medium-weight, multilayer non-woven material. It offers full protection against fluids under light pressure conditions.

KEY FEATURES

- AAMI PB70 level I protection.
- These gowns are for single-use.
- Side ties eliminate difficulties of tying behind the back.
- Neck ties provide a generous size for full coverage and flexibility
- They include functional cuffs constructed of a white knitted material.
- Provided latex-free.
- Provided non-sterile.

APPROVALS

8080 complies with the AAMI Protective Barriers classification system and the associated minimum requirements for the liquid barrier performance of protective apparel and drapes based on industry-accepted test methods.



SIZING

An appropriate size garment should be selected to allow sufficient movement for the task whilst maintaining a secure fit.

SIZE CHART

SIZE	SLEEVE LENGTH	CHEST WIDTH	BODY LENGTH
S / M	85 CM	73 CM	113 CM
L / XL	88 CM	76 CM	115 CM
XXL / 3XL	91 CM	82 CM	117 CM

MATERIALS

Fabric: PPSB nonwoven, 100% Polypropylene

Weight: 35-45 gsm

STORAGE AND DISPOSAL

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapours.
- Expected shelf life is three years from the date of manufacture when stored as stated.
- Replace garment if damaged, heavily contaminated or in accordance with local work practice or regulations.
- Handle and dispose of contaminated garments with care and in accordance with applicable regulations.

WARNING AND LIMITATIONS

Before use read and understand all user instructions and be sure that the product is suitable for the application and fitted correctly. Product must never be altered or modified. Final determination as to the suitability of these products for a particular situation is the employer's responsibility. This information is subject to revision at any time. If you have questions you may contact technical experts.

LIMITED USE



Single use



Do not wash



Do not bleach



Do not iron



Do not dry clean



Do not tumble dry



Flammable